miniNAV[®]

DE Gebrauchsanweisung | EN Instructions for use | ES Instrucciones de manejo | FR Mode d'emploi
| IT Istruzioni per l'uso | PT Instruções de utilização

 www.miethke.com

US This Instructions for Use is NOT intended for United States users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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TABLE OF CONTENTS

INDICATION	12
TECHNICAL DESCRIPTION	12
PHYSICAL BACKGROUND	12
FUNCTION OF THE VALVE	13
SELECTING THE APPROPRIATE SHUNT SYSTEM	14
POSSIBLE SHUNT COMPONENTS	14
TUBE SYSTEMS	15
IMPLANTATION	15
PREOPERATIVE VALVE TEST	16
RE-IMPLANTATION	16
SAFETY MEASURES AND CONTRAINDICATIONS	16
FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	16
POSTOPERATIVE VALVE TEST	16
ADVERSE REACTIONS	17
STERILISATION	17
RESTERILISATION	17
PRESSURE-FLOW CHARACTERISTICS	17
MEDICAL DEVICES CONSULTANTS	18
REQUIREMENTS OF THE MEDICAL DEVICE DIRECTIVE (DIRECTIVE 93/42/EEC)	18
DIMENSIONS	18

INDICATION

The *miniNAV* is used for draining cerebrospinal fluid from the ventricles into the peritoneum in hydrocephalus patients.

TECHNICAL DESCRIPTION

The *miniNAV* was developed to offer a small-sized valve without the obstruction-related problems that are known to arise in the treatment of hydrocephalus. The *miniNAV* is composed of a robust titanium casing (1) whose proximal end contains a ballcone valve. A spiral spring (2) maintains the opening pressure of the ball-cone valve and the sapphire ball (3) ensures the precise closure of the valve. The inlet connector (4) and the outlet connector (5) are also made of titanium.


Fig. 1: Schematic cross-section of the *miniNAV*

1. titanium housing
2. spiral spring
3. sapphire ball
4. inlet connector
5. outlet connector

PHYSICAL BACKGROUND

Horizontal position

When the body is in a lying position, the pressure at the valve is calculated by determining the difference between the intraventricular pressure and the pressure in the abdomen (see fig. 2). In a healthy person, the intraventricular pressure is positive when the body is in a lying position. To adjust this pressure by means of valve drainage, the appropriate pressure level must be chosen, taking the abdominal cavity pressure into account. The IVP is then calculated by determining the sum of the valve opening pressure and the abdominal cavity pressure (see fig. 2).

Vertical position

In a healthy person, the ventricular pressure is slightly negative when the body is in a standing position. To adjust this pressure by means of valve drainage, the valve opening pressure must be set much higher than would be required if the person was lying down. It is only in this way that the valve can compensate for the level of hydrostatic pressure minus the abdominal cavity pressure and the desirable, slightly negative, intraventricular pressure.


Fig. 2: Pressure ratios for lying and upright positions

IVP : Intraventricular pressure

P_{Vli} : Valve opening pressure in supine position (ball-cone valve only)

P_{Vst} : Valve opening pressure in standing position (ball-cone valve + gravitational valve)

P_B : Pressure in the abdominal cavity

P_{Hyd} : Hydrostatic pressure

Lying down: $IVP = P_{Vli} + P_B$

Standing: $IVP = P_{Hyd} - P_{Vst} - P_B$

FUNCTION OF THE VALVE

The operating principle of the miniNAV is illustrated in fig. 3. Fig. 3a shows the miniNAV in the horizontal position. The ball-cone valve is closed and drainage is prevented. If the patient's IVP increases and continues to rise, the spring pressure of the ball-cone unit will be overcome. The sealing ball will move away from the cone and a gap opens for fluid drainage.


Fig. 3: miniNAV

a) closed

b) open

SELECTING THE APPROPRIATE SHUNT SYSTEM

The *miniNAV* is available in four different pressure levels: 0, 5, 10 and 15 cmH₂O. The pressure setting should be chosen according to the clinical picture (normal-pressure hydrocephalus, hypertonic hydrocephalus).

In case the patient suffers from symptoms associated with overdrainage, or complications with overdrainage are expected, we recommend the implantation of the *SHUNTASSISTANT 2.0* in addition to the *miniNAV*. The *SHUNTASSISTANT 2.0* is a hydrostatic, supplementary valve specially designed for preventing problems with overdrainage. It is made by Christoph Miethke GmbH & Co. KG. The coding of the *miniNAV* can be identified according to the shape of the valve's housing. For example the *miniNAV* with an opening pressure of 5 cmH₂O has a concave proximal part (curved inwards) and a convex distal part (curved outwards) (fig. 4). Each *miniNAV* is calibrated in accordance with strict quality control standards. The following pressure levels are available:

Pressure rating (cmH ₂ O)	Coding
0	
5	
10	
15	


Fig. 4: X-ray image of the *miniNAV* (pressure rating 5 cmH₂O)

a. concave
b. convex

POSSIBLE SHUNT COMPONENTS

The *miniNAV* is available with different shunt accessories. These variations are comprised of a variety of components, which are introduced briefly below:

The *Burrhole Reservoir* is positioned in the cranial burr hole. It allows for measurement of the intraventricular pressure, injection of drugs and extraction of CSF. Its solid titanium base is highly puncture-resistant. All reservoirs are available with integrated catheters or connectors. A special *Burrhole Reservoir* is the *SPRUNG RESERVOIR*. An additional new feature of this reservoir is that CSF can be flushed towards the *miniNAV* because of a non-return valve in the bottom of the reservoir. By this mechanism, flow in the direction of the ventricular catheter is avoided during the pumping procedure. The opening pressure of the shunt system is not increased by the implantation of the *SPRUNG RESERVOIR*.

The *Prechamber* is positioned on the cranium. It allows for measurement of the intraventricular pressure, injection of drugs, extraction of CSF and palpatory ventricle inspection. Its solid titanium base is highly puncture-resistant. A puncture of the *Prechamber* or the *CONTROL RESERVOIR* should be performed as perpendicular to the reservoir surface as possible with a cannula of maximum diameter of 0.9 mm. 30 punctures are possible without any restrictions. A special *Prechamber* is the *CONTROL RESERVOIR*. As an additional new feature of this reservoir, CSF can be flushed towards the *miniNAV* because of a non-return valve in the

proximal inlet of the reservoir. By this mechanism, flow in the direction of the ventricular catheter is avoided during the pumping procedure. The opening pressure of the shunt system is not increased by the implantation of the *CONTROL RESERVOIR*.

WARNING

Frequent pumping can lead to overdrainage and thus to pressure conditions outside the normal physiological range. The patient should discuss the risks (involved) with their surgeon.

Tight tolerancing of the deflector ensures a good fit with the ventricular catheter. By adjusting the deflector (prior to implantation) the length of catheter penetrating into the skull can be optimised. The ventricular catheter is “deflected” at a right angle in the burr hole.

TUBE SYSTEMS

The *miniNAV* has been designed to ensure the optimal ventricular pressure. It is available as a shunt system or as individual valve units with or without an integrated distal catheter (internal diameter 1.2 mm, external diameter 2.5 mm). Individual valve units should be used with catheters of approx. 1.2 mm internal diameter and approx. 2.5 mm external diameter. The connector on the valve allows using catheters of 1.0 mm to 1.5 mm internal diameter. The external diameter of the catheter should be about double the internal diameter. Regardless, the catheters must be fixed carefully, with a ligature, to the valve connectors. It is essential that kinks in the catheter are avoided. The included catheters have virtually no effect on the pressure-flow characteristics.

IMPLANTATION

Positioning the ventricular catheter

Several surgical techniques are available for positioning the ventricular catheter. The necessary skin incision should be carried out, preferably, in the shape of a lobule pedicled towards the draining catheter or as a straight skin incision. To avoid CSF leakage, care should be taken that the dura opening is kept as small as possible after applying the burr hole. The ventricular catheter is stiffened by the introducing stylet supplied with the product.

The *miniNAV* is available in different shunt variations:

When using a *miniNAV* with a *Burrhole Reservoir* or *SPRUNG RESERVOIR*, the ventricular catheter is implanted first. Once the introducing stylet has been removed, the patency of the ventricular catheter can be tested by checking if CSF is dripping out. The catheter is shortened and the *Burrhole Reservoir* is connected, with the connection secured with a ligature. The skin incision should not be located directly above the reservoir. The *miniNAV* with *Prechamber* or *CONTROL RESERVOIR* comes with a deflector. This deflector is used for adjusting the position of deflection before implanting the ventricular catheter. The catheter is deflected; the *Prechamber* is put into place. The position of the ventricular catheter should be inspected again by postoperative CT or MR imaging.

Positioning the *miniNAV*

The *miniNAV* should be implanted in the head of the patient. The valve is marked with an arrow pointing to distal (downwards) to indicate the flow direction. Whether the label faces towards the skin or the brain is of no importance in terms of the valve's performance. Following subcutaneous tunneling, the catheter is either pushed from the burr hole, possibly through a reservoir, to the selected valve implantation site; or it is pushed through from the valve and connected to the reservoir, if there is any.

Positioning the peritoneal catheter

The access site for the peritoneal catheter is left to the surgeon's discretion. It can be applied e. g. para-umbilically in a horizontal direction or at the height of the epigastrium.

Likewise, various surgical techniques are available for positioning the peritoneal catheter. We recommend pulling through the peritoneal catheter, using a subcutaneous tunneling tool and perhaps with an auxiliary incision, from the shunt to the intended position of the catheter. The peritoneal catheter, which is usually attached securely to the *miniNAV*, has an open distal end, but no wall slits. Following the exposure of and the entry into the peritoneum by means of a trocar, the peritoneal catheter (shortened if necessary) is pushed forward into the open space in the abdominal cavity.

PREOPERATIVE VALVE TEST


Fig. 5: Patency test

a) Sterile, isotonic, physiological saline solution

The *miniNAV* can be filled by aspiration through a sterile, single-use syringe attached to the distal end of the catheter. The proximal end of the valve is immersed in a sterile, isotonic, physiological saline solution. The valve is patent if fluid can be extracted in this way (see fig. 5).

WARNING

Pressure admission through the single-use syringe should be avoided, both at the proximal and the distal end. Contaminations in the solution used for the test can impair the product's performance.

RE-IMPLANTATION

Products that have previously been implanted must not subsequently be reimplanted into the same or another patient.

SAFETY MEASURES AND CONTRAINDICATIONS

The patients must be carefully monitored after the implantation. Reddened skin and tension in the area of the drainage tissue could indicate infections at the shunt system. Symptoms such as headache, dizzy spells, mental confusion or vomiting are common occurrences in cases of shunt dysfunction. Such symptoms, as well as shunt system leakage, necessitate the immediate replacement of the shunt component responsible, or of the entire shunt system.

The implantation of medical devices is contraindicated if the patient has an infection or suspected infection (e. g. meningitis, ventriculitis, peritonitis, bacteraemia, septicaemia) in the region affected by the implantation.

FUNCTIONAL SAFETY AND COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

These medical devices are constructed in such a way as to ensure their precise and reliable operation over long periods of time. However, we cannot guarantee that these medical devices will require replacement for medical or technical reasons. These medical devices are able to resist positive and negative pressures up to 100 cmH₂O during and after implantation. These medical devices have to be stored in a clean and dry environment at all times.

Nuclear magnetic resonance (MRI) examinations up to a field strength of 3 Tesla or computed tomography (CT) examinations can be performed without risk or impairment to the valve function. The valve is MR Conditional. The supplied catheters are MR Safe. Reservoirs, deflectors and connectors are MR Conditional.

The conditions for MR compatibility of the products can be found on our website: <https://www.miethke.com/downloads/>

POSTOPERATIVE VALVE TEST

The *miniNAV* has been designed as a safe and reliable unit even without the implantation of a pumping device. However, the inclusion of a *Prechamber* or a *Burrhole Reservoir* allows the shunt system to be tested by flushing or pressure measurements.

ADVERSE REACTIONS

In the treatment of hydrocephalus with shunts, the following complications may arise (as described in the literature): infections, blockages caused by protein and/or blood in the cerebrospinal fluid, over-/underdrainage or in very rare cases noise development.

Violent shocks from the outside (accident, fall etc.) may put the integrity of the shunt system at risk.

STERILISATION

The products are sterilised with steam under closely monitored conditions. The expiry date is printed on the wrapping of each individual product. Products taken from a damaged wrapping must not be used under any circumstances.

RESTERILISATION

The functional safety and reliability of resterilized products cannot be guaranteed, therefore resterilisation is not recommended.

PRESSURE-FLOW CHARACTERISTICS

The diagrams below show the pressure-flow characteristics for the pressure ratings in which the *miniNAV* is available.





Fig. 6: The total opening pressure refers to a reference flow of 5 ml/h. When the flowrates reach 20 ml/h, the opening pressures are approximately 1 to 2 cmH₂O higher.

MEDICAL DEVICES CONSULTANTS

In compliance with the European Directive concerning medical devices (Directive 93/42/EEC), Christoph Miethke GmbH & Co. KG has nominated medical devices consultants as contacts for all product-related questions.

You can contact our medical devices consultants via:

Tel. +49 331 62083-0

info@miethke.com

REQUIREMENTS OF THE MEDICAL DEVICE DIRECTIVE (DIRECTIVE 93/42/EEC)

The Medical Device Directive requires comprehensive documentation of the fate of medical devices used in humans, especially for implants. The individual implant identification number should therefore be recorded in the patient's medical records and patient card to ensure complete traceability. Translations of these instructions for use into additional languages can be found on our website (<https://www.miethke.com/downloads>).

DIMENSIONS


Fig. 7: Schematic representation of the miniNAV with its external dimensions

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany