miniNAV[®]

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Ⓧ This Instructions for Use is NOT intended for United State users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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
CONTENT

INDICATION	13
TECHNICAL DESCRIPTION	13
PHYSICS BACKGROUND	14
FUNCTION OF THE VALVE	15
SELECTING THE APPROPRIATE SHUNTSYSTEM	15
POSSIBLE SHUNT COMPONENTS	15
TUBE SYSTEMS	16
SURGICAL PROCEDURE	16
PREOPERATIVE VALVE TEST	17
SAFETY MEASURES	17
RE-IMPLANTATION	17
COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	17
POSTOPERATIVE VALVE TEST	17
STERILISATION	17
RESTERILISATION	17
PRESSURE-FLOW CHARACTERISTICS	18
FUNCTIONAL SAFETY	18
SIDE EFFECTS	18
NOTE ON THE INSTRUCTIONS FOR USE	18
REQUIREMENTS OF THE MDD 93/42/EEC	18
MEDICAL PRODUCTS CONSULTANT	19
GENERAL INFORMATION	19
VARIATIONS	20

INDICATION

The *miniNAV* is used for draining cerebrospinal fluid from the ventricles into the peritoneum in hydrocephalus patients.

TECHNICAL DESCRIPTION

The *miniNAV* was developed to offer a small-sized valve without the obstruction related problems that are known to arise in the treatment of hydrocephalus.

The *miniNAV* is composed of a robust titanium casing (1) whose proximal end contains a ball-cone valve. A spiral spring (2) maintains the opening pressure of the ball-cone valve and the sapphire ball (3) ensures the precise closure of the valve.

The inlet connector (4) and the outlet connector (5) are also made of titanium.


Fig. 1: Schematic cross-section of the *miniNAV*

PHYSICS BACKGROUND

The intraventricular pressure is positive in a healthy human in a horizontal position. To maintain this pressure through shunt drainage, one has to choose the appropriate pressure range, taking into account the abdominal cavity pressure. The resulting IVP is the sum of the shunt opening pressure and the abdominal cavity pressure (fig. 6).

In a healthy human, the ventricular pressure in the vertical position becomes slightly negative. To maintain this pressure by means of shunt drainage, the shunt opening pressure has to be significantly higher so that the shunt can compensate for the hydrostatic pressure minus the sum of the abdominal cavity pressure and the slightly negative intraventricular pressure. Conventional shunts open immediately as soon as the patient stands up, which can lead to critical overdrainage.

- IVP Intraventricular pressure
- PVli Opening pressure (horizontal)
- PVst Opening pressure (vertical)
- PB Pressure in the abdominal cavity
- PHyd Hydrostatic pressure

horizontal: $IVP = PVli + PB$
 vertical: $IVP = PHyd - PVst - PB$



Fig. 2: Pressure conditions in horizontal and vertical body positions

FUNCTION OF THE VALVE


The operating principle of the *miniNAV* is illustrated in fig. 3.

Fig. 3a shows the *miniNAV* in the horizontal position. The ball-cone valve is closed and drainage is prevented.

If the patient's IVP increases and continues to rise, the spring pressure of the ball-cone unit will be overcome. The sealing ball will move away from the cone and a gap opens for fluid drainage.

a)

b)

Fig. 3: *miniNAV*

a) closed and b) open

SELECTING THE APPROPRIATE SHUNT-SYSTEM

The *miniNAV* is available in four different pressure levels (0, 5, 10 and 15 cmH₂O). The pressure setting should be chosen according to the clinical picture (normal-pressure hydrocephalus, hypertonic hydrocephalus).

In case the patient suffers from symptoms associated with overdrainage, or complications with overdrainage are expected, we recommend implantation of the *SHUNTASSISTANT* in addition to the *miniNAV*. The *SHUNTASSISTANT* is a hydrostatic, supplementary valve specially designed for preventing problems with overdrainage. It is made by Christoph Miethke GmbH & Co. KG.

The coding of the *miniNAV* can be identified according to the shape of the valve's housing. For example the *miniNAV* with an opening pressure of 5 cmH₂O has a concave proximal part (curved inwards) and a convex distal part (curved outwards). (Fig. 4)

Each *miniNAV* is calibrated in accordance with strict quality control standards. The following pressure levels are available:

pressure rating (cmH ₂ O)	Coding
0	
5	
10	
15	


Fig. 4: X-ray image of the *miniNAV* (pressure rating 5 cmH₂O)

POSSIBLE SHUNT COMPONENTS

The *miniNAV* is available with different shunt accessories. These variations are comprised of a variety of components, which are described briefly introduced below:

The *borehole reservoir* is positioned in the cranial borehole. It allows measurement of intraventricular pressure, injection of drugs and extraction of CSF. Its solid titanium base is highly puncture-resistant. All reservoirs are available with integrated catheters or connectors. A special borehole reservoir is the *SPRUNG RESERVOIR*. An additional new feature of this reservoir is that CSF can be flushed towards the *miniNAV* because of a non-return valve in the bottom of the reservoir. By this mechanism, flow in the direction of the ventricular catheter is avoided during the pumping procedure. The opening pressure of the shunt system is not increased by the implantation of the *SPRUNG RESERVOIR*.

The *prechamber* is positioned on the cranium. It allows measurement of intraventricular pressure, injection of drugs, extraction of CSF and palpatory ventricle inspection. Its solid titanium base is highly puncture-resistant. A puncture of the *prechamber* or the *CONTROL RESERVOIR* should be performed as perpendicular to the reservoir surface as possible with a cannula of maximum diameter 0,9mm. 30 punctures are possible without any restrictions. A special prechamber is the *CONTROL RESERVOIR*. As an additional new feature of this reservoir, CSF can be flushed towards the *miniNAV* because of a non-return valve in the proximal inlet of the reservoir. By this mechanism, flow in the direction of the ventricular catheter is avoided during the pumping procedure. The opening pressure of the shunt system is not increased by the implantation of the *CONTROL RESERVOIR*.

Warning note: Frequent pumping can lead to overdrainage and thus to pressure conditions outside the normal physiological range. The patient should discuss the risks (involved) with their surgeon.

Tight tolerancing of the deflector ensures a good fit with the ventricular catheter. By adjusting the deflector (prior to implantation) the length of catheter penetrating into the skull can be optimised. The ventricular catheter is “deflected” at a right angle in the borehole (see chapter “Variations”).

TUBE SYSTEMS

The *miniNAV* has been designed to ensure the optimal ventricular pressure. It is available as a shunt system or as individual valve units with or without an integrated distal catheter (internal diameter 1.2 mm, external diameter 2.5 mm). Individual valve units should be used with catheters of approx. 1.2 mm internal diameter and approx. 2.5 mm external diameter. The connector on the valve allows using catheters of 1.0 mm to 1.5 mm internal diameter. The external diameter of the catheter should be about double the internal diameter. Regardless, the catheters must be carefully fixed, with a ligature, to the valve connectors. It is essential that kinks in the catheter are avoided.

The included catheters have virtually no effect on the Pressure-flow characteristics.

SURGICAL PROCEDURE

Positioning the ventricular catheter

Several surgical techniques are available for positioning the ventricular catheter. The necessary skin incision should be carried out, preferably, in the shape of a lobule pedicled towards the draining catheter or as a straight skin incision. To avoid CSF leakage, care should be taken that the dura opening is kept as small as possible after applying the borehole. The ventricular catheter is stiffened by the introducing stylet supplied with the product.

The *miniNAV* is available in different shunt variations:

When using a *miniNAV SHUNTSYSTEM* with borehole reservoir or *SPRUNG RESERVOIR*, the ventricular catheter is implanted first. Once the introducing stylet has been removed, the patency of the ventricular catheter can be tested by checking if CSF is dripping out. The catheter is shortened and the borehole reservoir is connected, with the connection secured with a ligature. The skin incision should not be located directly above the reservoir.

The *miniNAV SHUNTSYSTEM* with prechamber or *CONTROL RESERVOIR* comes with a deflector. This deflector is used for adjusting the position of deflection before implantation of the ventricular catheter. The catheter is deflected; the prechamber is put into place. The position of the ventricular catheter should be inspected again by postoperative CT or MR imaging.

Positioning the *miniNAV*

The *miniNAV* should be implanted in the head of the patient.

The valve is marked with an arrow pointing to distal (downwards) to indicate the flow direction. Whether the label faces towards the skin or the brain is of no importance in terms of the valve’s performance.

Following subcutaneous tunneling, the catheter is either pushed from the borehole, possibly through a reservoir, to the selected valve implantation site; or it is pushed through from the valve and connected to the reservoir, if there is any.

Positioning the peritoneal catheter

The access site for the peritoneal catheter is left to the surgeon’s discretion. It can be applied e. g. para-umbilically in a horizontal direction or

transrectally at the height of the epigastrium. Likewise, various surgical techniques are available for positioning the peritoneal catheter. We recommend pulling through the peritoneal catheter, using a subcutaneous tunneling tool and perhaps with an auxiliary incision, from the shunt to the intended position of the catheter. The peritoneal catheter, which is usually securely attached to the *miniNAV*, has an open distal end, but no wall slits. Following the exposure of, and the entry into, the peritoneum by means of a trocar, the peritoneal catheter (shortened, if necessary) is pushed forward into the open space in the abdominal cavity.

PREOPERATIVE VALVE TEST


Fig. 5: Patency test

The *miniNAV* can be filled by aspiration through a sterile, single-use syringe attached to the distal end of the catheter. The proximal end of the valve is immersed in a sterile, physiological saline solution. The valve is patent if fluid can be extracted in this way (see Fig. 5).

Caution: Pressure admission through the single-use syringe should be avoided, both at the proximal and the distal end.

Contaminations in the solution used for the test can impair the product's performance.

SAFETY MEASURES

The patients must be carefully monitored after the implantation. Reddened skin and tension in the area of the drainage tissue could indicate infections at the shunt system. Symptoms such as headache, dizzy spells, mental confusion or vomiting are common occurrences in cases of shunt dysfunction. Such symptoms, as well as shunt system leakage, necessitate the imme-

diate replacement of the shunt component responsible, or of the entire shunt system

RE-IMPLANTATION

Under no circumstances should products that have had previously been implanted in a patient be subsequently reimplanted in another, because a successful decontamination of the device cannot be reached without functional degradation.

COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

MRI examinations with field strengths of up to 3.0 tesla and CT examinations can be carried without endangering or impairing the functionality of the shunt. The *miniNAV* is MR Conditional (ASTM-F2503-08). All components are visible via X-ray. The provided catheters are MRI Safe. Reservoirs, deflectors and connectors are MR Conditional.

POSTOPERATIVE VALVE TEST

The *miniNAV* has been designed as a safe and reliable unit even without the implantation of a pumping device. However, the inclusion of a prechamber or a borehole reservoir allows the shunt system to be tested by flushing or pressure measurements.

STERILISATION





The products are sterilized with steam under closely monitored conditions. Double wrapping in sterile bags ensures sterility for a period of five years. The expiry date is printed on the wrapping of each individual product. Products taken from a damaged wrapping must not be used under any circumstances.

RESTERILISATION

The functional safety and reliability of resterilized products cannot be guaranteed, therefore re-sterilisation is not recommended.

PRESSURE-FLOW CHARACTERISTICS

The diagrams below show the pressure-flow characteristics for the pressure ratings in which the miniNAV is available.

The total opening pressure refers to a reference flow of 5 ml/h. When the flowrates reach 20 ml/h, the opening pressures are approximately 1-2 cmH₂O higher.

FUNCTIONAL SAFETY

The valves have been designed for long-term reliable and precise operation. Still, the possibility that the shunt system will need to be replaced for technical or medical reasons cannot be excluded.

The valve and the valve system are able to resist positive and negative pressure up to 200 cmH₂O during and after implantation.

NOTE ON THE INSTRUCTIONS FOR USE

The descriptions and explanations given in this document are based on the clinical experience available to date. It is for the surgeon to decide if surgical procedures should be changed according to his or her experience and to surgical practice.

REQUIREMENTS OF THE MDD 93/42/EEC

The MDD calls for the comprehensive documentation of the whereabouts of medical products that are applied in human beings, especially the whereabouts of implants. For this reason, the individual identification numbers of any implanted valves are to be noted in patients' records, so that in the event of any inquiries, the implant can be traced without any difficulties. Each valve is outfitted with a sticker for this purpose.

MEDICAL PRODUCTS CONSULTANT

In compliance with the requirements of the European law MDD 93/42/EEC, Christoph Miethke GmbH&Co.KG names medical product consultants as the individuals to be addressed with all queries concerning the products:

Dipl.-Ing. Christoph Miethke
Dipl.-Ing. Roland Schulz

Christoph Miethke GmbH & Co. KG
Ulanenweg 2
D-14469 Potsdam
Phone: +49(0) 7000 6438453 or
Phone: +49(0) 331 620 83 0
Fax: +49(0) 331 620 83 40
e-mail: info@miethke.com


Please address any enquiries to:
AESCULAP AG
Am Aesculap Platz
D-78532 Tuttlingen
Phone: +49 (0) 7461 95-0
Fax: +49 (0) 7461 95-26 00
e-mail: information@aesculap.de

Service address in the US
AESCULAP Inc.
Attn. AESCULAP Technical Services
615 Lambert Pointe Road
Hazelwood, MO, 63042

AESCULAP Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420

Distributor in the US/ Contact in Canada
AESCULAP Inc.
3773 Corporate Parkway
Center Valley, PA 18034
Phone: +1-800-282-9000
www.aesculapusa.com

GENERAL INFORMATION

Manufacturer	Christoph Miethke GmbH & Co. KG
Product name	<i>miniNAV</i>
Intended use	Treatment of hydrocephalus
Intended for one-time use (disposable)	
Store in a clean, dry place	
Schematic representation of the shunt with its external dimensions:	
 <p style="text-align: center;">14,7 mm</p> <p style="text-align: right;">2,8 mm</p>	
scale: 1:1	

VARIATIONS

The *miniNAV* is available as a single valve or as a shunt system comprising various components.


miniNAV


miniNAV SHUNTSYSTEM

miniNAV SHUNTSYSTEM with SPRUNG RESERVOIR or with borehole reservoir
(adult and pediatric)

miniNAV SHUNTSYSTEM with CONTROL RESERVOIR or with prechamber
(adult and pediatric)

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany