paediGAV[®]

Ⓛ Gebrauchsanweisung | Ⓞ Instructions for use | Ⓧ Mode d'emploi
Ⓝ Instrucciones de manejo | ⓘ Istruzioni per l'uso

Ⓤ This Instructions for Use is NOT intended for United State users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000.

A paper copy will be provided to you upon request at no additional cost.

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CONTENT

INDICATION	14
TECHNICAL DESCRIPTION	14
PHYSICAL BACKGROUND	15
FUNCTION OF THE VALVE	16
SELECTING THE APPROPRIATE VALVE	16
READING THE PRESSURE SETTING FROM AN X-RAY IMAGE	16
POSSIBLE SHUNT COMPONENTS	17
TUBE SYSTEMS	17
SURGICAL PROCEDURE	18
PREOPERATIVE VALVE TEST	18
RE-IMPLANTATION	18
SAFETY MEASURES	19
COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	19
POSTOPERATIVE VALVE TEST	19
FUNCTIONAL SAFETY	19
ADVERSE REACTION	19
STERILIZATION	19
RESTERILIZATION	19
PRESSURE-FLOW CHARACTERISTICS	19
MEDICAL PRODUCTS CONSULTANT	21
NOTE ON THE INSTRUCTIONS FOR USE	21
REQUIREMENTS OF THE MDD 93/42/EEC	21
GENERAL INFORMATION	21
VARIATIONS	22

INDICATION

The *paediGAV* is intended for use in cases of pediatric hydrocephalus for the drainage of cerebrospinal fluid from the ventricles into the peritoneum.

TECHNICAL DESCRIPTION

The *paediGAV* is a position-sensitive valve for the treatment of hydrocephalus. It comprises a ball-cone unit and a gravitational unit. By this design, physiological drainage can be assured independent of the physical position of the patient.

Fig. 1 shows a schematic cross section of the *paediGAV*. It consists of a sturdy titanium casing, which houses, in its proximal section, a well-ried ball-cone unit. A spiral spring (2) maintains the opening pressure of the ball-cone unit. The gravitational unit in the distal part consists of a tantalum ball (3), which defines the opening pressure of this valve, and a sapphire ball (4), which ensures its precise closure.

As an option, a connector or a silicone catheter is installed at the distal end of the valve. Coding rings (1) allow the identification of the pressure setting in an x-ray image. The *paediGAV* regulates the intraventricular pressure (IVP) for the patient being in a horizontal or in a standing position.


Fig. 1: Schematic cross section of the *paediGAV*

PHYSICAL BACKGROUND

The intraventricular pressure is positive in a healthy human in a horizontal position. To adjust this pressure through shunt drainage, one has to choose the appropriate pressure range, taking into account the abdominal cavity pressure. The resulting IVP is the sum of the shunt opening pressure and the abdominal cavity pressure (fig.2)

The ventricular pressure in a healthy human in a vertical position becomes slightly negative. To maintain this pressure by means of shunt drainage, the shunt opening pressure has to be significantly higher so that the shunt can compensate the hydrostatic pressure minus the sum of the abdominal cavity pressure and the slightly negative intraventricular pressure. Conventional shunts open immediately as soon as the patient stands up, which can lead to critical overdrainage.

<i>IVP</i>	<i>Intraventricular pressure</i>
<i>PVli</i>	<i>Opening pressure (horizontal)</i> <i>(ball-cone unit only)</i>
<i>PVst</i>	<i>Opening pressure (vertical)</i> <i>(ball-cone unit+ gravitational unit)</i>
<i>PB</i>	<i>Pressure in the abdominal cavity</i>
<i>PHyd</i>	<i>Hydrostatic pressure</i>

horizontal: $IVP = PVli + PB$
vertical: $IVP = PHyd - PVst - PB$


Fig. 2: The pressure conditions for horizontal and standing positions

FUNCTION OF THE VALVE

The operating principle of the *paediGAV* is illustrated in fig. 3 and fig. 4.

Fig. 3 a shows the *paediGAV* in a horizontal position. The ball-cone unit is closed and the drainage is blocked. In fig. 3 b, the *paediGAV* is shown when it is open. The increased IVP of the patient overcomes the spring force, which otherwise would keep the ball-cone unit closed. The closing ball is pushed away from the cone seal and a gap opens for fluid drainage. The closing ball is always open in the horizontal position and does not present any resistance to the drainage flow.


Fig. 3: *paediGAV* in horizontal position
a) closed b) open

As soon as the patient raises to a vertical position, the gravitational unit closes and the drainage is interrupted (fig. 4 a). In this position, the sum of the IVP and the hydrostatic pressure surpasses the spring force of the ball-cone unit. Therefore this valve is open now. But only when the sum of the IVP and the hydrostatic pressure exceeds the gravitational force on the tantalum ball in the gravitational unit, the drainage path is open again (fig. 4 b).


Fig. 4: *paediGAV* in a vertical position
a) closed b) open

SELECTING THE APPROPRIATE VALVE

The *paediGAV* operates in different modes depending on the physical position of the patient. It features two different opening pressures, one for the patient in a horizontal position and one for a vertical position.

Percutaneous readjustments of the valve opening pressure become unnecessary through this design, as a higher opening pressure in the vertical position counteracts any unintentional overdrainage, while the lower opening pressure required for when the patient is lying down rules out any underdrainage.

Horizontal position:

For the horizontal position, two different opening pressure settings are available (4 and 9 cmH₂O). Selection of the appropriate pressure rating should be made in accordance with indication (age of the patient).

Vertical position:

The opening pressure for the vertical body position depends upon the size of the patient (hydrostatic parameter). If an infant is being treated, a low pressure rating should be selected. If, on the other hand, a larger child is being treated, a *paediGAV* with a higher pressure rating should be applied. (see recommendation of pressure settings at www.miethke.com)

READING THE PRESSURE SETTING FROM AN X-RAY IMAGE

Each *paediGAV* is calibrated under strict quality control. The following pressure combinations are available:

Opening pressure		Coding
horizontal	vertical	
4 cmH ₂ O	24 cmH ₂ O	
9 cmH ₂ O	24 cmH ₂ O	
9 cmH ₂ O	29 cmH ₂ O	
4 cmH ₂ O	14 cmH ₂ O	
4 cmH ₂ O	19 cmH ₂ O	
9 cmH ₂ O	19 cmH ₂ O	

Recommendation of pressure levels:

Opening pressure		Recommendation
horizontal	vertical	
4 cmH ₂ O	24 cmH ₂ O	children up to 6 months
9 cmH ₂ O	24 cmH ₂ O	children from 6 months to 5 years
9 cmH ₂ O	29 cmH ₂ O	children over 5 Years

The pressure settings selected can be checked postoperatively on x-ray images, on which the corresponding coding is visible.


Fig. 5: Radiograph of a paediGAV (9/24 cmH₂O)

POSSIBLE SHUNT COMPONENTS

The *paediGAV* is available with different shunt accessories. These variants consist of a variety of components, which are described briefly below:

The *borehole reservoir* is positioned in the cranial borehole. It allows measuring the intraventricular pressure, injecting drugs and extracting CSF. Its solid titanium base is highly puncture-resistant. All reservoirs are available with integrated catheters or connectors. A special *borehole reservoir* is the *SPRUNG RESERVOIR*. As additional new feature of this reservoir CSF can be flushed towards the valve because of a one-way valve in the bottom of the reservoir. By this mechanism a flow in the direction of the ventricular catheter is avoided during the pumping procedure. The opening pressure of the shunt system is not

increased by the implantation of the *SPRUNG RESERVOIR*.

The *prechamber* is positioned on the cranium. It allows measuring the intraventricular pressure, injecting drugs, extracting CSF and performing a palpatory ventricle inspection. Its solid titanium base is highly puncture-resistant. A puncture of the *prechamber* or the *CONTROL RESERVOIR* should be performed as perpendicular to the reservoir surface as possible with a cannula of max. 0,9 mm. 30 times of punctures are able without any restrictions. A special *prechamber* is the *CONTROL RESERVOIR*. As an additional new feature of this reservoir, CSF can be flushed towards the valve because of a one-way valve in the proximal inlet of the reservoir. By this mechanism a flow in the direction of the ventricular catheter is avoided during the pumping procedure. The opening pressure of the shunt system is not increased by the implantation of the *CONTROL RESERVOIR*.

Warning note: Frequent pumping can lead to overdrainage and thus to unphysiological pressure conditions. The patient should be informed about the risk.

Due to its tight fit on the *ventricular catheter*, the *deflector* allows choosing the length of catheter penetrating into the skull prior to implantation. The *ventricular catheter* is deflected at a right angle in the borehole

TUBE SYSTEMS

The *paediGAV* has been designed to ensure the optimal ventricular pressure. It is available as a shunt system or as individual valve units with or without an integrated distal catheter (internal diameter 1.2 mm, external diameter 2.5 mm). Individual valve units should be used with catheters of approx. 1.2 mm internal diameter and approx. 2.5 mm external diameter.

The connector on the valve allows using catheters of 1.0 mm to 1.5 mm internal diameter. The external diameter of the catheter should be about double the internal diameter. In any case, the catheters must be carefully fixed, with a ligature, to the valve connectors. Kinks in the catheter have to be avoided.

The provided catheters have virtually no effect on the Pressure-flow characteristics.

SURGICAL PROCEDURE

Positioning the ventricular catheter

Several surgical techniques are available for positioning the ventricular catheter. The necessary skin incision should be carried out, preferably, in the shape of a lobule pedicled towards the draining catheter or as a straight skin incision. To avoid CSF leakage, care should be taken that the dura opening is kept as small as possible after applying the borehole. The ventricular catheter is stiffened by the introducing stylet supplied with the product.

The *paediGAV* is available in different shunt variants: When using a *paediGAV SHUNTSYSTEM with borehole reservoir*, the ventricular catheter is implanted first. Once the introducing stylet has been removed, the patency of the ventricular catheter can be tested by checking if CSF is dripping out. The catheter is shortened and the *borehole reservoir* is connected, with the connection secured with a ligature. The skin incision should not be located directly above the reservoir.

The *paediGAV SHUNTSYSTEM with prechamber* comes with a *deflector*. This *deflector* is used for adjusting the position of deflection before implantation of the ventricular catheter. The catheter is deflected; the *prechamber* is put into place. The position of the ventricular catheter should be inspected again by postoperative CT or MR imaging.

Positioning the valve:

The *paediGAV* operates in different modes depending on the patient's position. Hence it is important that the valve is implanted parallel to the body axis. A suitable implant site is behind the ear. Following the skin incision and the tunneling through the skin, the catheter is pushed forward to the implantation site chosen for the valve. If necessary, the catheter is shortened and fastened at the *paediGAV* by means of a ligature, taking care that the valve does not lie directly under the skin incision. The valve is marked with an arrow pointing towards the feet of the patient (towards the distal end or downwards respectively).

Positioning the peritoneal catheter

The access site for the peritoneal catheter is left to the surgeon's discretion. It can be applied

e. g. para-umbilically in a horizontal direction or transectally at the height of the epigastrium. Likewise, various surgical techniques are available for positioning the peritoneal catheter. We recommend pulling through the peritoneal catheter, using a subcutaneous tunneling tool and perhaps with an auxiliary incision, from the shunt to intended position of the catheter. The peritoneal catheter, which is usually securely attached to the *paediGAV*, has an open distal end, but no wall slits. Following the exposure of, and the entry into, the peritoneum by means of a trocar, the peritoneal catheter (shortened, if necessary) is pushed forward into the open space in the abdominal cavity.

PREOPERATIVE VALVE TEST


Fig. 6: Patency test

The *paediGAV* can be filled most gently by aspiration through a sterile, single-use syringe attached to the distal end of the catheter. The proximal end of the valve is immersed in a sterile, physiological saline solution. The valve is patent if fluid can be extracted in this way (see fig. 6).

Caution: Pressure admission through the single-use syringe should be avoided, both at the proximal and the distal end.

Contaminations in the solution used for the test can impair the product's performance.

RE-IMPLANTATION

Under no circumstances should products that have had previously been implanted in a patient be subsequently reimplanted in another, because a successful decontamination of the device cannot be reached without functional degradation.

SAFETY MEASURES

The patients must be carefully monitored after the implantation. Reddened skin and tension in the area of the drainage tissue could indicate infections at the shunt system. Symptoms such as headache, dizzy spells, mental confusion or vomiting are common occurrences in cases of shunt dysfunction. Such symptoms, as well as shunt system leakage, necessitate the immediate replacement of the shunt component responsible, or of the entire shunt system

COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

MRI examinations with field strengths of up to 3.0 tesla and CT examinations can be carried without endangering or impairing the functionality of the shunt. The *paediGAV* is MR Conditional (ASTM F2503-13). All components are visible via X-ray. The provided catheters are MRI Safe. Reservoirs, *deflectors* and connectors are MR Conditional.

POSTOPERATIVE VALVE TEST

The *paediGAV* has been designed as a safe and reliable unit even without the provision of a pumping device. However, there are ways of testing the unit if a shunt system with a *prechamber* or a *borehole reservoir* is used. Valve tests can be carried out by flushing or pressure measurements.

FUNCTIONAL SAFETY

The valves have been designed for long-term reliable and precise operation. Still, it cannot be excluded that the shunt system needs to be replaced for technical or medical reasons. The valve and the valve system are able to resist positive and negative pressure up to 200 cmH₂O during and after implantation.

ADVERSE REACTION

In the treatment of hydrocephalus with shunts, the following complications may arise (as described in the literature): infections, blockages caused by protein and/or blood in the cerebrospinal fluid, over/under drainage or in very rare cases, noise development.

Due to violent shocks from the outside (accident, fall, etc.) the integrity of the shunt may be endangered.

STERILIZATION


The products are sterilized with steam under closely monitored conditions. The double wrapping in sterile bags ensures sterility for a period of five years. The expiry date is printed on the wrapping of each individual product. Products taken from a damaged wrapping must not be used under any circumstances.

RESTERILIZATION


The functional safety and reliability of resterilized products cannot be guaranteed, therefore re-sterilisation is not recommended.

PRESSURE-FLOW CHARACTERISTICS


In the following, the pressure-flow characteristics for the available pressure settings of the *paediGAV* are shown.


paediGAV 4/19 cmH₂O


paediGAV 9/24 cmH₂O


paediGAV 4/24 cmH₂O

paediGAV 9/29 cmH₂O

paediGAV 9/19 cmH₂O

The total opening pressure refers to a reference flow of 5 ml/h. When the flowrates reach 20 ml/h, the opening pressures are approximately 1-2 cmH₂O higher.

MEDICAL PRODUCTS CONSULTANT

In compliance with the requirements of the European law MDD 93/42/EEC, Christoph Miethke GmbH&Co. KG names medical product consultants as the individuals to be addressed with all queries concerning the products:

Dipl.-Ing. Christoph Miethke
Dipl.-Ing. Roland Schulz
Michaela Funk-Neubarth

Christoph Miethke GmbH & Co. KG
Ulanenweg 2
14469 Potsdam · Germany
Phone: +49(0) 7000 MIETHKE or
Phone: +49(0) 331 620 83 0
Fax: +49(0) 331 620 83 40
e-mail: info@miethke.com

Please address any enquiries to:
AESCULAP AG
Am Aesculap-Platz
78532 Tuttlingen · Germany
Phone: +49 (0) 7461 95 0
Fax: +49 (0) 7461 95 26 00
e-mail: information@aesculap.de

Service address in the US
AESCULAP Inc.
Attn. AESCULAP Technical Services
615 Lambert Pointe Road
Hazelwood, MO, 63042

AESCULAP Repair Hotline
Phone: +1 (800) 214-3392
Fax: +1 (314) 895-4420

Distributor in the US/ Contact in Canada
AESCULAP Inc.
3773 Corporate Parkway
Center Valley, PA 18034
Phone: +1-800-282-9000
www.aesculapusa.com


NOTE ON THE INSTRUCTIONS FOR USE

The descriptions and explanations given in this document are based on the clinical experience available to date. It is for the surgeon to decide if surgical procedures should be changed according to his or her experience and to surgical practice.

REQUIREMENTS OF THE MDD 93/42/EEC

The MDD calls for the comprehensive documentation of the whereabouts of medical products that are applied in human beings, especially the whereabouts of implants. For this reason, the individual identification numbers of any implanted valves are to be noted in patients' records, so that in the event of any inquiries, the implant can be traced without any difficulties. Each valve is outfitted with a sticker for this purpose.

GENERAL INFORMATION

Manufacturer	Christoph Miethke GmbH & Co. KG
Product name	<i>paediGAV</i>
Intended use	Treatment of Hydrocephalus
Intended for one-time use (disposable)	
Store in a clean, dry place	
Schematic representation of the valve with its external dimensions:	
	

VARIATIONS

The *paediGAV* is available as a single valve or as a shuntsystem comprising various components.

paediGAV

paediGAV SHUNTSYSTEM

paediGAV SHUNTSYSTEM with pediatric borehole reservoir

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany