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Recommended pressure levels for GAV 2.0® and GAV 2.0® LP

The GAV 2.0 comprises a differential pressure unit and a gravitational unit. The opening pressure of the GAV 2.0 in the horizontal position is determined by the differential pressure unit. The opening pressure, which works in the vertical position is the sum of both: differential pressure unit and gravitational unit. With the help of the Miethke App it is easy to understand the technical functionality of the GAV 2.0 in different body positions.


Please download the Miethke App to see the operating principle of our gravitational valves.

Example // GAV 2.0:
5 / 25 cmH₂O

Horizontal Position:
Differential pressure unit with 5 cmH₂O

Vertical Position:
Differential pressure unit and Gravitational unit act together:
5 / 25 cmH₂O:
5 cmH₂O + 20 cmH₂O = 25 cmH₂O

Differential pressure unit

Gravitational unit

Values apply for mobile patients. In case of less mobile patients or a high BMI the gravitational unit should be set lower than it is recommended here.

This is a non-binding recommendation.

The physician decides in each case individually.

	Children from 6 months	Children from 3 years	Adults	Adults from 65 years
Newborn	10 / 25 *	10 / 30 *	5 / 30 *	5 / 25 *
			5 / 25 * < 1,60 m 5 / 35 * > 1,80 m	5 / 20 * < 1,60 m 5 / 30 * > 1,80 m

* recommended pressure level in cmH₂O

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