## *DUALSWITCH VALVE (DSV)®*

ⓓ Gebrauchsanweisung | Ⓒ Instructions for use | Ⓕ Mode d'emploi  
Ⓔ Instrucciones de manejo | ⓘ Istruzioni per l'uso

Ⓤ This Instructions for Use is NOT intended for United State users. Please discard.  
The Instructions for Use for United States users can be obtained by visiting our website at [www.aesculapusa.com](http://www.aesculapusa.com) and clicking the "Products" menu. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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
## CONTENT

INDICATION	16
TECHNICAL DESCRIPTION	16
FUNCTION OF THE DSV	16
PHYSICAL BACKGROUND	18
SELECTING THE APPROPRIATE DSV SHUNTSYSTEM	19
READING THE PRESSURE SETTING FROM AN X-RAY IMAGE	19
POSSIBLE SHUNT COMPONENTS	20
TUBE SYSTEMS	20
SURGICAL PROCEDURE	20
TESTING THE VALVE PATENCY	22
PRESSURE-FLOW-CHARACTERISTICS	22
RE-IMPLANTATIONS	23
SAFETY MEASURES	23
COMPATIBILITY WITH DIAGNOSTIC PROCEDURES	23
POSTOPERATIVE VALVE TESTS	24
FUNCTIONAL SAFETY	24
ADVERSE REACTION	24
STERILISATION	24
RESTERILISATION	24
REQUIREMENTS OF THE MDD 93/42/EEC	24
NOTE ON THE INSTRUCTIONS FOR USE	24
MEDICAL PRODUCTS CONSULTANT	25
GENERAL INFORMATION	25
VARIANTS	26

**INDICATION**

The *DUALSWITCH VALVE* is used for fluid drainage from the ventricles into the peritoneum in cases of hydrocephalus.

**TECHNICAL DESCRIPTION**

The *DUALSWITCH VALVE* was developed with the aim to avoid known problems (blocked drainage system, overdrainage, influence of the subcutaneous pressure) in the treatment of hydrocephalus. The *DUALSWITCH VALVE* features two parallel chambers, one for the patient vertical and one for the horizontal position. Fig. 1 shows a schematic cross section of the *DUALSWITCH VALVE*. The robust titanium casing (1) contains two titanium discs (3) connected to the casing by two silicone membranes (2). Two yoke springs (5), (6) of different strength act on these discs and determine the opening pressure with their respective preload force. The springs press the pressure discs against a ball (4) integrated into the casing. In this way, the drill holes in the discs are closed and the valve is closed. The stronger one of the two

springs (5) controls the opening and closing for the high-pressure side (7); the weaker spring (6) determines the opening and closing of the low-pressure side (8). A secondary tantalum ball (9) on the lowpressure side acts as the switch between the two pressure levels. Inlet port (10) and outlet port (11) are also made of titanium.

**FUNCTION OF THE *DUALSWITCH VALVE***

The *DUALSWITCH VALVE* regulates the intraventricular pressure in different ways depending on the position of the patient (vertical or horizontal).

The operating principle of the *DUALSWITCH VALVE* is illustrated in figures 2 and 3. Fig. 2a shows the closed valve in the lying-down position. Both membranes are pressed onto the casing ball; the valve blocks any fluid drainage. Fig. 2b shows the valve open (lying-down position). The IVP of the patient has risen, the spring pressure on the low-pressure side is overcome, and the membrane is pressed up from the ball, thereby opening a gap for fluid drainage. The pressure cannot rise any further.


Fig. 1: Schematic cross section of the *DUALSWITCH VALVE*


Fig. 2: Valve states in the lying-down position: a) closed, b) open

As soon as the patient rises, the tantalum ball closes the drainage channel on the low-pressure side. As long as the sum of the IVP and the hydrostatic pressure does not exceed the closing force of the stronger spring, the valve remains closed, see fig. 3a. Only if the pressure rises further, the high-pressure side of the valve opens and thereby prevents any further pressure rise, see fig. 3b.


Fig. 3: Valve states in the vertical position: a) closed, b) open

The *DUALSWITCH VALVE* has large active switching surfaces, which ensures precise operation. The valve opening force, F, which is brought about by ventricular pressure, is governed by the level of ventricular pressure, p, and the surface, A, on which the ventricular pressure acts. In conventional valves this active surface is very small. Consequently, the valve opening forces are small as well.

Thus, countervailing forces (sticking caused by fluid residues, e.g. blood cells) can severely impair the valve's functionality.

$$p = \frac{F}{A} = \frac{F_{interference}}{A} = \frac{F}{A} = \frac{F}{A}$$

Conventional Valve
*DUALSWITCH-Valve*

To minimize the effects of such interfering forces, but still be able to use the same valve pressure rating, the opening force needs to be increased. In this way, the sticking forces will be very small compared to the opening force so that their influence on valve performance is reduced considerably. Assuming the ventricular pressure is unchanged, the opening force can be increased by providing for a larger active surface in the valve design. In the *DUALSWITCH VALVE*, the fluid is distributed over a surface significantly larger than in conventional valves. Therefore the valve is much less susceptible to conglutinations of the membrane through proteins (fibrin)-rich cerebrospinal fluid.

To realize the advantage of strong opening forces, the valve is about 28 x 5mm in size. The *DUALSWITCH VALVE* is placed in the thoracic region, like cardiac pacemakers, which have been implanted there successfully for many years.


Fig. 4: Fluid and pressure distribution in the closed valve

**PHYSICAL BACKGROUND**

In the lying-down position, the intraventricular pressure is positive in a healthy human. To adjust this pressure using valve drainage, only the pressure in the abdominal cavity has to be taken into account.

Fig. 5 Calculation of the intraventricular pressure for the horizontal and standing body position. In the standing body position, the ventricular

pressure should be slightly negative. To regulate ventricular pressure using shunt valve drainage, valve opening pressure in standing position must be set considerably higher. Only in this way can the shunt valve compensate for the hydrostatic pressure reduced by the pressure in the abdominal cavity. Conventional valves open instantly, as soon as the patient rises to a standing position, and a critical overdrainage may occur.


Fig. 5: Pressure situation at the shunt for horizontal and vertical position

**SELECTING THE APPROPRIATE**

*DSV SHUNTSYSTEM*

The mode of operation of the *DUAL SWITCH-VALVE* depends on the position of the patient. It features two different opening pressures. One opening pressure is effective when the patient is in a horizontal position; the other is for the patient standing upright. Percutaneous readjustment of the valve characteristics is unnecessary, since the higher opening pressure in the upright position counteracts any unintentional overdrainage, while the lower opening pressure in the horizontal position prevents underdrainage.

**Horizontal position:**

For the horizontal position, four different opening pressure ratings are available (5, 10, 13 and 16 cmH<sub>2</sub>O). The recommended standard pressure rating is 10 cmH<sub>2</sub>O (valve pressure at a flow of 20 ml/h). The pressure rating 5 cmH<sub>2</sub>O is available for normal pressure hydrocephalus patients.

**Vertical position:**

The opening pressure for the standing position depends on the body height of the patient (hydrostatics). The high-pressure side of the valve is chosen in such a way that, with the patient standing upright, a ventricular pressure of at least -5 cmH<sub>2</sub>O is maintained under all circumstances. The procedure is the following:


- 1) The distance between the third ventricle (at the level of the Foramen Monroi) and the diaphragm of the patient is measured.
- 2) 5 cm are subtracted from the measured distance.
- 3) The valve is chosen whose pressure rating for the high-pressure side exceeds the measured distance minus 5 cm by the smallest amount. The available opening pressure ratings are 30, 40 and 50 cmH<sub>2</sub>O. The ventricular pressure in the patient will then be kept between -5 cmH<sub>2</sub>O and +5 cmH<sub>2</sub>O at all times.

**READING THE PRESSURE SETTING FROM AN X-RAY IMAGE**


Strict quality control is applied in the calibration of each individual *DUALSWITCH VALVE*. The following pressure level combinations, which can also be read from an X-ray visible coding, are available:

The code for the low-pressure side is always situated close to the ball switch (white dot). The code for the high-pressure side is always further away from the ball switch, no matter on which side of the ball switch the code is visible.


horizontal 5 cmH<sub>2</sub>O:


Standard-valve: horizontal 10 cmH<sub>2</sub>O:

horizontal 13 cmH<sub>2</sub>O:

horizontal 16 cmH<sub>2</sub>O:


Fig. 6 Possible codes (top picture) and an X-ray image of the DUALSWITCH VALVE (pressure ratings 13/30)

In case of questions or problems concerning the valve functionality, please contact the medical products advisors at Christoph Miethke GmbH & CO. KG.

## POSSIBLE SHUNT COMPONENTS

The DSV is available with different shunt accessories. These variants are comprised of a variety of components, which are described briefly introduced below:

The *Borehole Reservoir* is positioned in the cranial borehole. It allows measuring the intraventricular pressure, injecting drugs and extracting CSF. Its solid titanium base is highly puncture-resistant. All reservoirs are available with integrated catheters or connectors. A special *Borehole Reservoir* is the *SPRUNG RESERVOIR*. As additional new feature of this reservoir CSF can be flushed towards the valve because of a one-way valve in the bottom of the reservoir. By this mechanism a flow in the direction of the ventricular catheter is avoided during the pumping procedure. The opening pressure of the shunt system is not increased by the implantation of the *SPRUNG RESERVOIR*.

The prechamber is positioned on the cranium. It allows measuring the intraventricular pressure, injecting drugs, extracting CSF and performing a palpatory ventricle inspection. Its solid titanium base is highly puncture-resistant. A puncture of the *Prechamber* or the *CONTROL RESERVOIR* should be performed as perpendicular to the reservoir surface as possible with a cannula of max. 0,9 mm. 30 times of punc-

tures are able without any restrictions. A special *Prechamber* is the *CONTROL RESERVOIR*. As an additional new feature of this reservoir, CSF can be flushed towards the valve because of a one-way valve in the proximal inlet of the reservoir. By this mechanism a flow in the direction of the ventricular catheter is avoided during the pumping procedure. The opening pressure of the shunt system is not increased by the implantation of the *CONTROL RESERVOIR*.

**Warning note: Frequent pumping can lead to overdrainage and thus to unphysiological pressure conditions. The patient should be informed about the risk.**

## TUBE SYSTEMS

The DSV has been designed to ensure the optimal ventricular pressure. It is available as a shunt system or as individual valve units with or without an integrated distal catheter (internal diameter 1.2 mm, external diameter 2.5 mm). Individual valve units should be used with catheters of approx. 1.2 mm internal diameter and approx. 2.5 mm external diameter. The connector on the valve allows using catheters of 1.0 mm to 1.5 mm internal diameter. The external diameter of the catheter should be about double the internal diameter. This rule does not apply where the valve is used for lumboperitoneal drainage. In this case, a lumbar catheter can be connected to the valve or the reservoir, respectively. In any case, the catheters must be carefully fixed, with a ligature, to the valve connectors. Kinks in the catheter have to be avoided. The provided catheters have virtually no effect on the pressure-flow characteristics.

## SURGICAL PROCEDURE

### Positioning the ventricular catheter

Several surgical techniques are available for positioning the ventricular catheter. The necessary skin incision should be carried out, preferably in the shape of a lobule pedicled towards the

draining catheter or as a straight skin incision. To avoid CSF leakage, care should be taken that the dura opening is kept as small as possible after applying the borehole. The ventricular catheter is stiffened by the introducing stylet supplied with the product.

The DSV is available in different shunt variants: When using a *DSV SHUNTSYSTEM* with *Borehole Reservoir* or *SPRUNG RESERVOIR*, the ventricular catheter is implanted first. Once the introducing stylet has been removed, the patency of the ventricular catheter can be tested by checking if CSF is dripping out. The catheter is shortened and the *Borehole Reservoir* is connected, with the connection secured with a ligature. The skin incision should not be located directly above the reservoir.

The *DSV SHUNTSYSTEM* with prechamber or *CONTROL RESERVOIR* comes with a deflector. This deflector is used for adjusting the position of deflection before implantation of the ventricular catheter. The catheter is deflected; the prechamber is put into place. The position of the ventricular catheter should be inspected again by postoperative CT or MR imaging.

### Positioning a lumbar catheter

Presently, a lumbar catheter is not part of the product range of Christoph Miethke GmbH & CO. KG. The relevant manufacturer's instructions must be observed when using a lumbar catheter.

### Positioning the DSV

The mode of operation of the *DUALSWITCH VALVE* depends on the body position of the patient. Therefore, care must be taken to ensure that the connector axis is implanted parallel to the body axis, as far as possible.

The implantation is performed subcutaneously, either supraclavicular or just infraclavicular. The valve is marked with an arrow pointing distally (downward) to indicate the flow direction. The valve is implanted flat in a subcutaneous pocket; the orientation of the labeling (towards the skin or towards the ribs) does not affect the valve functionality.

Following the incision and the tunneling under the skin, the catheter is either pushed from the borehole to the chosen valve implantation site or the other way round, depending on which system is used. For lumboperitoneal drainage, we recommend choosing the arcus costalis as the implant position.

**Warning note: The adjustable unit must be placed over a hard boney surface and should not be implanted within an area, which makes finding and feeling the valve more difficult (e. g. under a scar).**

**Frequent pumping can lead to overdrainage and thus to unphysiological pressure conditions. The patient should be informed about the risk.**

### Positioning the peritoneal catheter

The access site for the peritoneal catheter is left to the surgeon's discretion. It can be applied e. g. para-umbilically in a horizontal direction or transrectally at the height of the epigastrium. Likewise, various surgical techniques are available for positioning the peritoneal catheter. We recommend pulling through the peritoneal catheter, using a subcutaneous tunneling tool and perhaps with an auxiliary incision, from the shunt to intended position of the catheter.

The peritoneal catheter, which is usually securely attached to the DSV, has an open distal end, but no wall slits. Following the exposure of, and the entry into, the peritoneum by means of a trocar, the peritoneal catheter (shortened, if necessary) is pushed forward into the open space in the abdominal cavity.

**TESTING THE VALVE PATENCY**

The DSV can be filled most gently by aspiration through a sterile, single-use syringe attached to the distal end of the catheter. The proximal end of the valve is immersed in a sterile, physiological saline solution. The valve is patent if fluid can be extracted in this way (see fig. 7).














Fig. 7: Patency test

**Caution:** Pressure admission through the single-use syringe should be avoided, both at the proximal and the distal end. Contaminations in the solution used for the test can impair the product's performance.

**PRESSURE-FLOW CHARACTERISTICS**

The diagrams below show the pressure-flow characteristics for the pressure ratings in which the *DUALSWITCH VALVE* is available.

The total opening pressure refers to a reference flow of 5 ml/h. When the flowrates reach 20 ml/h, the opening pressures are approximately 1-2 cmH2O higher.

**RE-IMPLANTATION**

Under no circumstances should products that have had previously been implanted in a patient be subsequently reimplanted in another, because a successful decontamination of the device cannot be reached without functional degradation.

**SAFETY MEASURES**

The patients must be carefully monitored after the implantation. Reddened skin and tension in the area of the drainage tissue could indicate infections at the shunt system. Symptoms such as headache, dizzy spells, mental confusion or vomiting are common occurrences in cases of shunt dysfunction. Such symptoms, as well as shunt system leakage, necessitate the immediate replacement of the shunt component responsible, or of the entire shunt system.

**COMPATIBILITY WITH DIAGNOSTIC PROCEDURES**

MRI examinations with field strengths of up to 3.0 tesla and CT examinations can be carried out without endangering or impairing the functionality of the shunt. The DSV is MR conditional (ASTM-F2503-13). All components are visible via X-ray. The provided catheters are MRI safe. Reservoirs, deflectors and connectors are MR conditional. Warning note: When using a magnetic field and simultaneous pressing on the valve an adjustment of the valve cannot be excluded.

**POSTOPERATIVE VALVE TEST**

The *DUALSWITCH VALVE* was designed as a safe and reliable unit operating without a pump or test equipment, so as to prevent unwanted interference by unauthorized or unqualified persons. However, there are ways of testing if the *Prechamber* or the *Borehole Reservoir* is used. Apart from that, the functioning of the tantalum ball switch can be checked by means of X-ray imaging of the valve with the patient standing upright and lying down. The tantalum ball, which switches between the two valve levels, appears as a near-white circle in the X-ray picture. In the picture for the horizontal position, a dark shadow must be visible directly below the ball, in the direction of the fluid flow; this shade is not visible in the picture for the upright position. A shadow being visible independent of the patient's position (oblong bar, see fig. 8b) means that the tantalum ball does not close when the patient stands up. In this case, the valve operates like a conventional differential pressure valve. If no shadow is visible, independent of the patient's position, the valve operates in high-pressure mode even when the patient is lying down, which means that the valve must be replaced.


Fig. 8: *DUALSWITCH VALVE* in a) vertical and b) horizontal position

**FUNCTIONAL SAFETY**

The valves have been designed for long-term reliable and precise operation. Still, the possibility that the shunt system needs to be replaced for technical or medical reasons, cannot be excluded. The valve and the valve system are able to resist positive and negative pressure up to 200 cmH<sub>2</sub>O during and after implantation.

**Warning note for carriers of pacemakers: Due to the implantation of a DSV the function of a pacemaker can be affected.**

**ADVERSE REACTION**

In the treatment of hydrocephalus with shunts, the following complications may arise (as described in the literature): infections, blockages caused by protein and/or blood in the cerebrospinal fluid, over/under drainage or in very rare cases, noise development.

Due to violent shocks from the outside (accident, fall, etc.) the integrity of the shunt may be endangered.

**STERILISATION**

The products are sterilized with steam under closely monitored conditions. The double wrapping in sterile bags ensures sterility for a period of five years. The expiry date is printed on the wrapping of each individual product. Products taken from a damaged wrapping must not be used under any circumstances.

**RESTERILISATION**

The functional safety and reliability of resterilized products cannot be guaranteed, therefore resterilisation is not recommended.

**REQUIREMENTS OF THE MDD 93/42/EEC**

The MDD calls for the comprehensive documentation of the whereabouts of medical products that are applied in human beings, especially the whereabouts of implants. For this reason, the individual identification numbers of any implanted valves are to be noted in patients' records, so that in the event of any inquiries, the implant can be traced without any difficulties. Each valve is outfitted with a sticker for this purpose.

**NOTE ON THE INSTRUCTIONS FOR USE**

The descriptions and instructions given above are based on clinical experience available to date. The shunt systems have been implanted according to the instructions given in this document, in clinical trials. It is the responsibility of the individual surgeon to change procedures according to her/his experience and surgical practice. He or she also does so at their own risk. However, the valve should always be implanted in the thorax area (or, in case of L-P drainage, on the arcus costalis), and not in the region of the head. Implantation in the thorax area ensures that the valve switches to the appropriate pressure mode when the patient changes his or her position.

**MEDICAL PRODUCTS CONSULTANT**


In compliance with the requirements of the European law MDD 93/42/EEC, Christoph Miethke GmbH&Co.KG names medical product consultants as the individuals to be addressed with all queries concerning the products:

Dipl.-Ing. Christoph Miethke  
 Dipl.-Ing. Roland Schulz  
 Christoph Miethke GmbH & Co. KG  
 Ulanenweg 2  
 14469 Potsdam · Germany  
 Phone: +49(0) 7000 MIETHKE or  
 Phone: +49(0) 331 620 83 0  
 Fax: +49(0) 331 620 83 40  
 e-mail: info@miethke.com

Please address any enquiries to:  
 AESCULAP AG  
 Am Aesculap-Platz  
 78532 Tuttlingen · Germany  
 Phone: +49 (0) 7461 95 0  
 Fax: +49 (0) 7461 95 26 00  
 e-mail: information@aesculap.de

**GENERAL INFORMATION**

<b>Manufacturer</b>	<b>Christoph Miethke GmbH &amp; Co. KG</b>
Product name	<i>DUALSWITCH VENTIL</i>
Intended use	Treatment of Hydrocephalus
Intended for one-time use (disposable)	
Store in a clean, dry place	


## VARIANTS

The DSV is available as a single valve or as a shunt system comprising various components.


### DUALSWITCH VALVE


### DUALSWITCH SHUNTSYSTEM (without reservoir)


### DUALSWITCH SHUNTSYSTEM with Borhole Reservoir or SPRUNG RESERVOIR (adult and pediatric)

### DUALSWITCH SHUNTSYSTEM with Prechamber or CONTROL RESERVOIR (adult and pediatric)

### DUALSWITCH SHUNTSYSTEM for L-P drainage

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany