## *miniNAV*<sup>®</sup>

Ⓓ Patientenhandbuch | ⒸⒺ ⒺⒶ Patient Manual | Ⓕ ⒻⒹ Manuel du patient  
ⒺⒹ Manual para el paciente | ⒺⒹⒹ Manuale per il paziente

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---

## THE COMPANY

Christoph Miethke GmbH & Co. KG is a company based in Berlin-Brandenburg region that develops, manufactures and markets innovative neurosurgical implants for the treatment of hydrocephalus. In the course of our work we established successful partnerships with numerous hospitals worldwide. The purpose of this booklet is to provide you and your family with some understanding of the treatment of hydrocephalus. The successful treatment of this condition has only been possible since the 1950s. In a dramatic race against time to save the life of his son, Casey, who suffered from hydrocephalus, a technician named John D. Holter developed, in only a few weeks, a novel silicone valve. Despite the fact that, since its first implantation in March 1956, this valve has proven to be clinically effective and a giant step in the treatment of this condition, there are many patients today who experience considerable problems with hydrocephalus valve systems. Christoph Miethke utilised the knowledge gained in 50 years of valve treatment and developed a new generation of highprecision valves made of the metal titanium.

## BASIC ANATOMICAL FACTS

The human brain (Fig. 1) is surrounded by a special substance known as cerebrospinal fluid (CSF). Cerebrospinal fluid is produced in several chambers, so-called ventricles, that are found within the brain. The channels, by which the ventricles are interconnected, constitute a complex drainage system. The fluid in the brain circulates through these ventricles and eventually flows into the venous blood. The function of this fluid is to protect the brain from mechanical damage. The CSF also regulates the internal brain pressure (intracranial pressure, ICP), keeps the brain tissue moist and transports the products of metabolism.

Fig. 1 Anatomy of the skull  
(see inner cover page)

- 1) Cranium
- 2) Brain
- 3) Cerebrospinal fluid
- 4) Lateral ventricle
- 5) Third ventricle
- 6) Fourth ventricle

## CLINICAL PICTURE

In healthy people, a balance exists between the production and resorption of cerebrospinal fluid. In infants, approximately 100 ml of this fluid are produced, whereas in small children about 250 ml is produced and in adults 500 ml. If the amount of fluid produced exceeds the amount resorbed, the ventricles expand, leading to the disorder known as hydrocephalus. Hydrocephalus refers to a condition whereby the volume of the "water" (hydro) in the head (cephalus) continually increases. This condition is often observed at birth (congenital hydrocephalus). The disorder can also develop later in life, e.g., as the result of inflammation, a severe head injury, a cancerous growth or scalp inflammation. Such cases are referred to as acquired hydrocephalus. A distinction is made between hydrocephalus occlusus (obstructive hydrocephalus) and hydrocephalus communicans (communicating hydrocephalus). In hydrocephalus occlusus, the links between the ventricles of the brain are severed in such a way that they are unable to "communicate" with each other. Cases in which the channels between the ventricles are unobstructed but resorption of cerebrospinal fluid is impaired are defined as hydrocephalus communicans.


Fig. 2 Size of the ventricle  
a) normal, b) hydrocephalus

## CLINICAL SYMPTOMS

In infants, the cranial bones have not yet grown together solidly. The increasing volume of cerebrospinal fluid causes the head to increase in circumference while, at the same time, brain tissue disintegrates. From the age of about two, the

hardened skull prevents any growth of the head's circumference. In this case, the increase in fluid volume leads to a massive pressure increase, resulting in the expansion of the brain ventricles and the compression of the brain itself. The consequence for infants and adults can be irreversible brain damage. Symptoms (depending on the severity of the disorder) include nausea, headache, vomiting, impaired coordination, drowsiness and, in the end, unconsciousness.

## DIAGNOSIS OF THE CONDITION

Doctors have a variety of ways at their disposal to diagnose hydrocephalus. The ventricle size is measured through imaging procedures (e.g. computerized tomography, ultrasound or NMR-tomography).

### Computerized tomography (CT)

This quick diagnostic procedure produces X-ray images of different layers of the head.

### Nuclear Magnetic Resonance (NMR) tomography

This electromagnetic imaging process produces images of very fine layers of the head. It is also known as NMR, MRT, or MRI scanning.

### Ultrasound

This procedure, in which the interior of the head is examined through the open fontanel, can only be applied to small children.

Another way of diagnosing hydrocephalus is through pressure measurements showing an increased brain pressure. The circulation of cerebrospinal fluid is investigated through examinations with contrast agents.

## TREATMENT MODALITIES

For all the efforts to find therapeutic alternatives to valve implantation (e. g. through pharmaceutical treatment or, most recently, by minimally invasive surgery), there is currently no alternative, in most cases, to the implantation of a drainage system, referred to as a shunt.

- 1 Right atrium
- 2 Cardiac catheter (atrial catheter)
- 3 *miniNAV*
- 4 Reservoir
- 5 Ventricular Catheter
- 6 Ventricles of the brain
- 7 Abdominal catheter (peritoneal catheter)
- 8 Abdominal cavity


Fig. 3 Drainage systems for hydrocephalus patients  
a) ventriculo-atrial, b) ventriculo-peritoneal

## COMPLICATIONS ARISING FROM THERAPY

A number of different complications can arise from shunt therapy for hydrocephalus, including infections or undesirable obstruction of the drainage system. The obstruction can result from increased protein concentrations in the cerebrospinal fluid. The *miniNAV* prevents such complications from occurring. This is achieved through a sufficiently large sectional area within the valve.

An additional complication is an undesirably high level of drainage of cerebrospinal fluid while the patient is in standing position. This can be avoided through the implantation, in addition to the *miniNAV*, of a gravitational valve - the *SHUNTASSISTANT*.

The *SHUNTASSISTANT* is a hydrostatic supplementary valve that was developed specifically to prevent overdrainage. It is a product of Christoph Miethke GmbH & Co.KG.

## POSTOPERATIVE BEHAVIOR

As a rule, the everyday activities of patients with shunt implants are not restricted. However, patients should consult their attending physician before major physical exertion (e. g. hard physical work, strenuous sports). Hydrocephalus patients who experience headache, dizziness, unnatural gait or similar symptoms should consult a physician without delay. Apart from that, we recommend medical check-ups at regular intervals. The patient should avoid knocks or pressure on the valve and catheters. The valve has been designed to be resistant against magnetic fields.

## VALVE MECHANISM

The *miniNAV* is outfitted with a conical ball valve that provides optimal valve closure and accurate pressure regulation. The sturdy titanium housing ensures that the influence of external factors such as subcutaneous pressure are negligible. This guarantees a high level of functional reliability. The unique design and its compact construction makes the *miniNAV* particularly suitable for both newborns and adults. Like other valves, the *miniNAV* is implanted on the cranium in such a way that the patient's activities will not be restricted by the presence of the implant. The valve can also easily and safely be implanted in other regions of the body and under the skin. Your attending physician will determine which implantation site is best for you.


Fig. 5 *miniNAV*, Scale 1:1

Fig. 4 Schematic cross section of the *miniNAV* (see inner cover page)

- 1 inlet connector
  - 2 titanium casing
  - 3 sapphire ball
  - 4 coil spring
  - 5 outlet connector
- a ball-cone-valve

**Warning note:** The shunt system may comprise a reservoir that can be pumped. Frequent pumping can lead to overdrainage and thus pressure conditions outside of the normal physiological range. As such pumping should only be carried out by the physician.

## PATIENT ID

Every *miniNAV* comes with a patient ID. The ID card is filled out by the attending physician and contains important information pertaining to follow-up examinations.

## A BRIEF PATIENT GLOSSARY

### Anatomy

A guide to the structure of body components

### Arachnoid

Connective tissue in the brain that lies between the dura mater and the pia mater

### Catheter

Tube

### Cerebrospinal fluid

Watery spinal fluid in the brain

### Communicating vessels

Vessels that are connected by a channel

### Computed tomography (CT)

Imaging technique whereby „slices“ of the body are recorded with an X-ray scanner

### Drainage

Drainage of accumulated fluid

### Dura mater

The hardest component of the meninx

### Fluid component

Cerebrospinal fluid component

### Fontanel

Connective tissue opening in a young infant's skull that later ossifies

### Implant

Substance that is placed in the human body to replace a particular function for a limited period of time or for the rest of the patient's life

### Lumbar puncture

Puncture of the spinal channel at the lower spine

### Lumbo-peritoneal drainage

Drainage of cerebrospinal fluid from the ventricle of the brain, by way of the region of the lumbar vertebrae in the abdominal cavity

### Meninges

Membrane found in the brain and spine

### Meningitis

Inflammation of the meninx

### Minimally invasive

Minimally infiltrating

### Overdrainage

Undesirable outward flow of cerebrospinal fluid

### Peritoneum

Membrane that covers the pelvic and abdominal cavities

### Pia mater

Component of the soft meninx containing blood vessels

### Pia arachnoid

Soft component of the meninx that is divided into the arachnoidea and pia mater

### Puncture

Insertion of a hollow needle or a trocar into a vessel for the purpose of removing fluid

### Resorption

Suctioning or removal of material through skin, mucosa or tissue

### Shunt

A passage between two channels – here a catheter drainage system with an integrated valve

### Spinal column

Element of the central nervous system located within the vertebral channel

### Subcutaneous pressure

Pressure beneath the skin

### Subdural hematoma

An accumulation of blood between the brain and cranium

### Ventricle of the brain

Intracranial space containing cerebrospinal fluid

### Ventricular peritoneal drainage

Drainage of cerebrospinal fluid from the ventricle of the brain directly into the abdominal cavity (abdominal catheter)

---

**FOLLOW-UP EXAMINATIONS**

A follow-up examination must be carried out in all cases.

Date	Treatment

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany