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

CHRISTOPH MIETHKE GMBH & CO. KG

Adjustable DP (differential pressure) –unit	Gravitational unit		
Standard (children and NPH-patients)	5 cmH ₂ O	Children up to 5 years	20 cmH ₂ O
Defensive (e.g. patients with extremely wide ventricles and highly elevated ICP or aqueductal stenosis)	10 cmH ₂ O	Children (over 5 years) & adults (up to 60 years)	25 cmH ₂ O
Special (e.g. patients with pseudotumor cerebri)	15 cmH ₂ O	Adults over 60 years	20 cmH ₂ O

These recommendations are based on common patient treatments, but can vary depending on the individual patient's condition.

Overweight

The peritoneal pressure inhibits drainage. Therefore the use of the gravitational unit should be considered for overweight patients as a function of body-mass-index (BMI):

- 25 - 29 BMI - 5 cmH₂O below standard recommendation
- 30 - 34 BMI - 10 cmH₂O below standard recommendation
- 35 - 39 BMI - 15 cmH₂O below standard recommendation
- > 40 BMI only apply the adjustable DP-unit (no gravitational unit necessary)

Mobility

Standard pressure levels are suitable for active people. Bedridden patients should not be treated with a gravitational unit.

Height

The hydrostatic suction effect normally depends on the height. Therefore we recommend the following corrections for the gravitational unit:

- < 1,60 m height: - 5 cmH₂O
- > 1,80 m height: + 5 cmH₂O

Control

The adjustable DP-unit is adjusted to 5 cmH₂O by the manufacturer. Please control the pressure setting before use.

Body position

Only the adjustable DP-unit determines the pressure level in the horizontal position. The analogous pressure level for the vertical body position results from the addition of both - adjustable DP-unit and gravitational unit.

Adjustment

Changing the adjustable DP-unit also influences the opening pressure level in the vertical body position.

■ CHRISTOPH MIETHKE GMBH & CO. KG