proSA[®] Tools

(DE) Gebrauchsanweisung | (EN) Instructions for use | (FR) Mode d'emploi | (ES) Instrucciones de manejo
(PT) Instruções de utilização | (IT) Istruzioni per l'uso | (NL) Gebruiksaanwijzing

 www.miethke.com

(US) This Instructions for Use is NOT intended for United States users. Please discard.

The Instructions for Use for United States users can be obtained by visiting our website at www.aesculapusa.com. If you wish to obtain a paper copy of the Instructions for Use, you may request one by contacting your local Aesculap representative or Aesculap's customer service at 1-800-282-9000. A paper copy will be provided to you upon request at no additional cost.

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Tel. +49 331 62083-0
info@miethke.com

TABLE OF CONTENTS

0.00 PREFACE AND IMPORTANT INFORMATION	22
1.00 INFORMATION ON HANDLING THESE INSTRUCTIONS FOR USE	22
1.01 Explanation of the warnings	22
1.02 Display conventions	22
1.03 Other accompanying documents and additional information	22
1.04 Feedback on the instructions for use	23
1.05 Copyright, disclaimer, warranty and other information	23
2.00 DESCRIPTION OF THE <i>proSA Tools</i>	23
2.01 Medical purpose	23
2.02 Clinical use	23
2.03 Indications	23
2.04 Contraindications	23
2.05 Intended patient groups	23
2.06 Intended users	23
2.07 Intended usage environment	23
2.08 Technical description	24
2.09 System components	25
3.00 PROPERTIES OF THE <i>proSA Tools</i>	26
3.01 Product description	26
3.02 Important safety information	27
3.03 Transport and storage	27
3.04 Use of the product	27
3.05 Disposal	33
3.06 Troubleshooting	33
3.07 Technical information	35
3.08 Symbols used for labelling	36
4.00 MEDICAL DEVICE CONSULTANTS	36

0.00 PREFACE AND
IMPORTANT INFORMATION

Preface

Thank you for purchasing *proSA Tools*. Please contact us if you have any questions about the contents of these instructions for use or the use of the product.

Your team at Christoph Miethke GmbH & Co. KG.

Relevance of the instructions for use

**WARNING**

Improper handling and use of this product can cause risks and damage. Therefore, please read and closely follow these instructions for use. Always keep them to hand. Follow the safety instructions to avoid personal injury or material damage.

Scope

The *proSA Tools* are comprised of the following components:

- ▶ *proSA Verification Tool*
- ▶ *proSA Masterdisc*
- ▶ *proSA Verification Compass*
- ▶ *proSA Adjustment Tool*
- ▶ *proSA Adjustment Disc*

Additional options:

- ▶ *proSA Checkmate*

1.00 INFORMATION ON HANDLING THESE
INSTRUCTIONS FOR USE

1.01 EXPLANATION OF THE WARNINGS

**DANGER**

Indicates an immediate risk. If not avoided, it will lead to death or serious injuries.

**WARNING**

Indicates a possible risk. If not avoided, it could lead to death or serious injuries.

**CAUTION**

Indicates a possible risk. If not avoided, it may lead to slight or minor injuries.

**NOTE**

Indicates a potentially harmful situation. If not avoided, the product or items in its vicinity may be damaged.

The symbols for risk, warning and caution are yellow warning triangles with a black edge and a black exclamation mark.

1.02 DISPLAY CONVENTIONS

Display	Description
<i>Italics</i>	Indicates product names

1.03 OTHER ACCOMPANYING
DOCUMENTS AND ADDITIONAL
INFORMATION

These instructions for use as well as translations into additional languages can be found on our website:

<https://www.miethke.com/downloads/>

If you still require help despite carefully reading the instructions for use and the additional information, please contact us or your authorised distributor.

1.04 FEEDBACK ON THE INSTRUCTIONS FOR USE

Your opinion is important to us. Please let us know if you have any requests and criticisms about these instructions for use. We will analyse your feedback and take it into account for the next version of the instructions for use where appropriate.

1.05 COPYRIGHT, DISCLAIMER, WARRANTY AND OTHER INFORMATION

Christoph Miethke GmbH & Co. KG guarantees a faultless product that is free of material and manufacturing defects upon delivery.

No liability, guarantee or warranty for safety and functionality can be assumed if the product is modified in any way other than described in this document, if it is combined with products by another manufacturer or if it is used in any way other than for the intended purpose and the intended use.

Christoph Miethke GmbH & Co. KG points out that the reference to its trademark rights applies solely to jurisdictions in which it has trademark rights.

2.00 DESCRIPTION OF THE *proSA Tools*

2.01 MEDICAL PURPOSE

The *proSA Tools* are used to set and check the adjustable pressure levels of the *proSA* valve.

2.02 CLINICAL USE

Treatment options

- ▶ Identifying the pressure level of the *valve's Gravitational Unit proSA*
- ▶ Checking the pressure level without diagnostic procedures such as X-rays
- ▶ Setting the pressure level without invasive procedures

2.03 INDICATIONS

The *proSA Tools* are used for the following indications:

- ▶ Treatment of hydrocephalus

2.04 CONTRAINDICATIONS

The following contraindications apply to the *proSA Tools* :

- ▶ Intolerance to materials of the instruments

2.05 INTENDED PATIENT GROUPS

Patients who on account of their clinical picture are treated with a CSF-draining shunt system with an adjustable *proSA* valve.

2.06 INTENDED USERS

In order to avoid risks due to false diagnoses, incorrect treatments and delays, the product must only be used by users with the following qualifications:

- ▶ Medical specialists, e.g. neurosurgeons
- ▶ Knowledge about the function and proper use of the product
- ▶ Successful participation in product training

2.07 INTENDED USAGE ENVIRONMENT

Professional healthcare facilities

- ▶ Healthcare facilities, such as hospitals, doctors' surgeries or rehabilitation facilities

2.08 TECHNICAL DESCRIPTION

proSA Verification Tool

The *Verification Tool* (Fig. 1) is used to read the set valve opening pressure. The *Verification Tool* is centred over the valve on the skin. It contains two magnets. As soon as the display (1) is released by pressing the push-button (2), the magnets of the pin line up with the magnets in the valve and the opening pressure is shown. The marking on the *Verification Tool* (3) must face the proximal catheter (leading to the ventricle) when being placed on the valve. The arrow on the bottom of the cap indicates the direction of flow of the shunt.

Fig. 1: proSA Verification Tool

proSA Masterdisc

Prior to reading the valve opening pressure, the *Verification Tool* can be checked by being placed onto the *proSA Masterdisc* (Fig. 2). The positions 0, 10, 20, 30 and 40 cmH₂O are indicated on the *proSA Masterdisc* (Fig. 2). If the display window is lined up with the value on the *Masterdisc* when placing the *Verification Tool*, the two values must match.

Example: The *Verification Tool* is placed on the *Masterdisc* so that the display window points to the value 10 on the *Masterdisc*. In this case, the *Verification Tool* should also show the value 10.

Fig. 2: proSA Masterdisc

proSA Verification Compass

In addition to the *Verification Tool*, there is a further option for measuring the set opening pressure. The *proSA Verification Compass* (Fig. 3) is placed onto the skin over the implanted valve and moved in small circles until the inner disc is firmly aligned. The opening pressure corresponds to the value displayed in the direction of the *Ventricular Catheter*.

NOTE

Small air bubbles in the *Verification Compass* chamber do not impair its function.

Fig. 3: proSA Verification Compass

proSA Adjustment Tool

The *proSA Adjustment Tool* (Fig. 4) is used to set the valve opening pressure. The desired pressure level is selected by turning the adjusting wheel (1); the pressure level appears on the display window (2). The *Adjustment Tool* is centred over the valve. Pressing the push-button (3) presses the adjustment tip (4) out of the base, unlocks the rotor in the valve and sets the valve opening pressure. The marking (5) on the *Adjustment Tool* must face the proximal catheter (leading towards the ventricle) when being placed onto the valve. The arrow on the bottom of the cap indicates the direction of flow.


Fig. 4: *proSA Adjustment Tool*

proSA Adjustment Disc

The pressure level can also be changed using the *proSA Adjustment Disc* (Fig. 5). The *Adjustment Disc* is centred over the valve. When doing so, the desired pressure level must point in the direction of the *Ventricular Catheter*. Slightly pressing the *Adjustment Disc* onto the valve releases the brake and changes the pressure level of the adjustable Gravitational Unit. The *proSA Adjustment Disc* has a red housing and is only available in one size. The scale shows the available pressure levels from 0 to 40 cmH₂O.

Fig. 5: *proSA Adjustment Disc*
Colour: red

2.09 SYSTEM COMPONENTS

The *proSA Tools* can be used to set the opening pressures of the following valves:

Product name	
proSA	

3.00 PROPERTIES OF THE *proSA Tools*

3.01 PRODUCT DESCRIPTION


3.01.01 APPROVED COMPONENTS

The *proSA Tools* are comprised of the following components:

- ▶ *proSA Adjustment Tool* (1),
- ▶ *proSA Verification Tool* (2),
- ▶ *proSA Masterdisc* (3),
- ▶ *proSA Verification Compass* (4) and
- ▶ *proSA Adjustment Disc* (5).

The *proSA Tools* can also be supplemented with the

- ▶ *proSA Checkmate* (6).

	Item no.	Figure
1	FV790T	
2	FV791T	
3	FV793T	
4	FV794T	
5	FV795T	
6	FV796T	

3.01.02 SCOPE OF DELIVERY

Box content (item no. FV792T)	Number
<i>proSA Adjustment Tool</i>	1
<i>proSA Verification Tool</i>	1
<i>proSA Masterdisc</i>	1
<i>proSA Verification Compass</i>	1
<i>proSA Adjustment Disc</i>	1
<i>proSA Instruments</i> Instructions for use	1

Box content (item no. FV796T)	Number
<i>proSA Checkmate</i>	1
<i>proSA Instruments</i> Instructions for use	1

3.01.03 STERILITY

With the exception of the *proSA Checkmate*, the *proSA Tools* cannot be sterilised, only disinfected.

WARNING

If the packaging or the product is damaged, the *proSA Tools* must not be used.

3.01.04 REPEATED USE AND RE-STERILISATION

The *proSA Checkmate* must be sterilised in double sterile packaging in an autoclave (steam sterilisation, fractionated vacuum method) at 134 °C and 5 minutes retention time.

3.01.05 OPERATING CONDITIONS

The *proSA Tools* are to be used at room temperature in healthcare facilities, such as hospitals, doctors' surgeries and rehabilitation facilities.

3.01.06 PRODUCT LIFE

These medical devices are designed to operate accurately and reliably over long periods of time. The expected product life of the *proSA Tools* is three years after first use, provided that the product is exposed to normal usage conditions and is properly maintained (see chapter "3.04 Use of the product").

The instruments may be used beyond this time if they are in perfect working condition. How-

ever, we cannot guarantee that these medical devices will require replacement for medical or technical reasons.

3.01.07 PRODUCT CONFORMITY

The product meets current regulatory requirements.

Translations of these instructions for use into additional languages can be found on our website:

<https://www.miethke.com/downloads/>

3.02 IMPORTANT SAFETY INFORMATION

3.02.01 SAFETY INSTRUCTIONS

Important! Read all safety instructions carefully before using the product. Follow the safety instructions in order to avoid injuries and life-threatening situations.

WARNING

- ▶ If the sterile packaging of the Checkmate is damaged, the product must not be used.
- ▶ The Checkmate must not be used after the expiry date and must be re-sterilised as described in chapter 3.01.04.
- ▶ Due to the risk of injury resulting from incorrect use of the *proSA Tools*, the instructions for use must be read carefully and understood before the product is used for the first time.

WARNING

- ▶ Prior to use, it is essential to check that the *proSA Tools* are complete and intact.

3.02.02 COMPLICATIONS, SIDE EFFECTS AND RESIDUAL RISKS

Various side effects and complications can occur when using the *proSA Tools*:

- ▶ Allergic reaction / intolerance to the materials of the instruments

Immediately consult a physician if the patient suffers from skin rashes and tightness, severe headaches, dizzy spells or similar.

3.02.03 REPORTING OBLIGATION

All serious incidents (damage, injuries, infections, etc.) occurring in relation to the product must be reported to the manufacturer and the responsible state authority.

3.02.04 PATIENT INFORMATION

The attending physician is responsible for informing the patient and/or his/her proxy in advance. The patient is to be informed about warnings, precautions, contraindications, precautionary measures to be taken as well as restrictions on use in relation to the product.

3.03 TRANSPORT AND STORAGE

The medical devices must always be transported and stored in a clean and dry place.

Storage conditions

Temperature range for storage	0 °C to +40 °C
-------------------------------	----------------

3.04 USE OF THE PRODUCT

3.04.01 INTRODUCTION

The *proSA Tools* must only be used by trained specialists.

The *proSA Tools* can be used to determine, change and monitor the pressure level of the adjustable Gravitational Unit of the *proSA* valve.

3.04.02 SAFETY AND WARNING NOTICES

CAUTION

- ▶ The *proSA Adjustment Tool* emits a magnetic field. Metallic objects and magnetic storage media should be placed at a sufficient safety distance.
- ▶ Because of the magnets inside the *proSA Tools*, the instruments must not be used in the vicinity of cardiac pacemakers. Furthermore, there is risk of damage to MRI instruments in their vicinity. For that reason, the *proSA Tools* must not be used in such locations!

CAUTION

- It is absolutely essential that only the *proSA Tools* are used to determine, change and monitor the opening pressure of the gravitational unit of the *proSA*.

- If the *proSA* is combined with an adjustable differential pressure unit (*proGAV*), *proGAV Tools* must be used to adjust and check the adjustable pressure level of the *proGAV* valve.

3.04.03 REQUIRED MATERIALS

To adjust the adjustable Gravitational Unit of the *proSA* after implantation, the first requirement is for the *proSA Verification Tool* / the *proSA Verification Compass* in order to read the set valve opening pressure, followed by the *proSA Adjustment Tool*/the *proGAV Adjustment Disc* to adjust the valve opening pressure.

The *proSA Masterdisc* can be used to check the Verification Tool.

To adjust the adjustable Gravitational Unit of the *proSA* before or during valve implantation, the resterilisable *proSA Checkmate* can be used.

3.04.04 PREPARING FOR USE

Checking the sterile packaging

Immediately before using the *proSA Checkmate*, the sterile packaging must be visually inspected in order to check the integrity of the sterile barrier system. The products should only be removed from the packaging immediately prior to use.

Checking the integrity of the *proSA Tools*

Prior to use, the instruments must be checked for integrity. To this end, all instruments must be visually inspected.

3.04.05 USING THE INSTRUMENTS

Adjusting the adjustable Gravitational Unit of the *proSA* requires the following steps:

1. Localisation

The adjustable Gravitational Unit is located under the skin (Fig. 6).


Fig. 6: Localisation of the adjustable Gravitational Unit

2. Test procedure

NOTE

The *proSA Verification Tool* and the *proSA Verification Compass* react sensitively to external magnetic fields. In order to rule out unwanted interactions, the *proSA Adjustment Tool* should not be placed in the immediate vicinity of the *proSA Verification Tool* or the *proSA Verification Compass* when determining the opening pressure. We recommend a distance of approx. 30 cm.

2a. Testing with the *proSA Verification Tool*

WARNING

The *Verification Tool* must be centred over the valve if possible, otherwise the determined opening pressure may be incorrect!

The *Verification Tool* is centred over the valve. The marking on the *Verification Tool* must face the proximal catheter (leading to the ventricle) (Fig. 7).


Fig. 7: Placing the Verification Tool

The release button on the *Verification Tool* is pressed and the set pressure level is then read. The arrow on the bottom of the *Verification Tool* indicates the direction of flow of the cerebrospinal fluid.


Fig. 8: Pressure level measurement

a) correct

b) wrong

2b. Testing with the proSA Verification Compass

In addition to the *Verification Tool*, there is a further option for measuring the set opening pressure. The *proSA Verification Compass* (Fig. 3) is placed onto the skin over the implanted valve and moved in small circles until the inner disc is firmly aligned. The opening pressure corresponds to the value displayed in the direction of the *Ventricular Catheter*.

NOTE

Small air bubbles in the *Verification Compass* chamber do not impair its function.

3. Adjustment process

WARNING

When adjusting the adjustable Gravitational Unit, care must be taken that the opening pressure is changed by a maximum of 16 cmH₂O per adjustment process.

Example: The opening pressure for the adjustable Gravitational Unit on the *proSA* is to be changed from 6 to 36 cmH₂O. If there is only one adjustment process, the rotor will rotate anti-clockwise (short route) and stop at the setting 0 cmH₂O (Fig. 9a).

The correct method is an adjustment in two stages: Initial adjustment from 6 cmH₂O to 22 cmH₂O and then from 22 cmH₂O to 36 cmH₂O. The rotor rotates clockwise (Fig. 9b).


Fig. 9: Rotor rotation during adjustment

a) wrong direction

b) correct direction

NOTE

- During preoperative adjustment through the sterile packaging, the *Adjustment Tool* is only placed on the valve. Do not press the push-button. Excessive pressure with the *Adjustment Tool* can damage the valve housing and impair its function.
- If the implantation site is unfavourably chosen or the skin above the valve is too thick, the adjustable Gravitational Unit may no longer be able to be adjusted. The valve then works with unchangeable pressure levels. The patient should be made aware of this risk.
- Swelling of the skin may make adjustment difficult for a few days after the operation.

3a. Adjustment with the proSA Adjustment Tool

WARNING

When pressing the release button, make sure that the *Adjustment Tool* is not removed from the valve.

NOTE

For sensitive patients, the use of local anaesthetic should be considered during the adjustment process (e.g. with a plaster), provided there is no contraindication.

Turning the adjusting wheel until the desired opening pressure is shown on the display sets the *Adjustment Tool* to the desired opening pressure (Fig. 10).


Fig. 10: Setting the desired pressure level on the *Adjustment Tool*

The *Adjustment Tool* is centred over the valve. The marking (a)/the display window must face the proximal catheter (leading to the ventricle) (b) (Fig. 11).


Fig. 11: Placing the *Adjustment Tool*

Once the *Adjustment Tool* is centred on the valve, press the push-button and the adjustment tip will come out of the base of the *Adjustment Tool*. This triggers the mechanical brake in the valve and sets the desired opening pressure (Fig. 12).


Fig. 12: Setting the pressure level

3b. Adjustment with the proSA Adjustment Disc

To use the *Adjustment Disc* for adjustment, it must be centred on the valve so that the desired pressure level (b) points in the direction of the inlet connector or the *Ventricular Catheter* (c) (Fig. 13).


Fig. 13: Adjustment with the *Adjustment Disc*

Press on the valve with the *Adjustment Disc* (Fig. 14) without changing its position. Then remove the *Adjustment Disc* and check the pressure level with the *Verification Tool*.


Fig. 14: Placing and pressing the Adjustment Disc onto the adjustable Gravitational Unit

3c. Adjustment with the proSA Checkmate $\text{C}\text{€}$ 0297

The *Checkmate* is supplied sterile and can be re-sterilised. This makes it possible to carry out a pressure level change and the check during valve implantation directly on the adjustable Gravitational Unit.

To determine the pressure level, the *Checkmate* is centred over the adjustable Gravitational Unit. The *Checkmate* automatically aligns itself over the valve. The pressure level can be read off towards the proximal catheter (leading to the ventricle).

To adjust the pressure level, the *Checkmate* is centred over the adjustable Gravitational Unit. When doing so, the required pressure level must point towards the proximal catheter (leading to the ventricle). By slightly pressing the *Checkmate* onto the valve, the adjustable Gravitational Unit brake is released and the pressure level is set.

Fig. 15: proSA Checkmate

Colour: gold

Pressure levels: 0 to 40 cmH₂O

4. Checking after adjustment

A check is carried out after adjusting the valve opening pressure. To do this, proceed as in Point 2. If the measured pressure does not match the required pressure level, the adjustment process should be repeated. To do this, start again at Point 3.

NOTE

If the valve setting cannot be checked conclusively using the *Verification Tool* or the *Verification Compass*, we recommend carrying out a check using an imaging method.

3.04.06 CLEANING RECOMMENDATION FOR NON-STERILISABLE proSA Tools

NOTE

The *proSA Tools* are made from thermolabile components that are sensitive to heat and moisture as well as chemically reactive. Do not immerse the *proSA Tools* in cleaning solutions or allow liquids to enter the housing, as moisture, corrosion and impurities could impair proper functioning.

Surface contamination on the *proSA Tools* should always be removed immediately after use by wiping with alcohol-based cleaning agents (minimum 75% alcohol).

The contact time should be at least 60 seconds and should be adapted to the degree of contamination. The instruments should then be wiped with a dry cloth.

The following cleaning methods are not suitable for processing the *proSA Tools* (except for the *Checkmate*): Irradiation, ultrasound, sterilisation, mechanical processing, immersion in cleaning fluids.

**3.04.07 CLEANING AND DISINFECTING
THE *proSA Checkmate***

Avoid damaging the product due to unsuitable cleaning agents/disinfectants and/or too high temperatures!

- ▶ Use cleaning agents/disinfectants approved for surgical steels in accordance with the manufacturer's instructions.
- ▶ Observe specifications regarding concentration, temperature and contact time.
- ▶ Do not exceed the maximum permitted cleaning temperature of 55°C.

Perform ultrasonic cleaning:

- ▶ as an effective mechanical support for manual cleaning/disinfection.
- ▶ for precleaning products with dried residues prior to mechanical cleaning/disinfection.
- ▶ as integral mechanical support during mechanical cleaning/disinfection.
- ▶ for the subsequent cleaning of products with residues still present after mechanical cleaning/disinfection.
- ▶ Provided the instruments can be fastened securely and safely for cleaning in machines or on storage aids, clean and disinfect the instruments mechanically.

Manual cleaning/disinfection

- ▶ Check visible surfaces for residues after manual cleaning/disinfection.
- ▶ If necessary, repeat the cleaning process.

Mechanical cleaning/disinfection

- ▶ Place product onto suitable screen basket.

Mechanical cleaning/disinfection with manual precleaning

Manual precleaning with ultrasound

Phase	I	II
Step	Disinfecting ultrasonic cleaning	Intermediate rinsing
T (°C/°F)	RT (cold)	RT (cold)
t (min)	15	1
Conc. (%)	2	-
Water quality	D-W	D-W
Chemical specifications	B. Braun Stabimed®; free of aldehydes, phenols and QAV; pH = 9	

D-W: drinking water; RT: room temperature

Mechanical alkaline cleaning and thermal disinfection

Device type: Single chamber cleaning/disinfection device without ultrasound

- Place product onto suitable screen basket (avoid rinse shadows)

Phase I

- Clean product in the ultrasonic cleaning bath (frequency 35 kHz). Ensure that all accessible surfaces are wetted.

Phase II

- Fully rinse/flush product (all accessible surfaces) under running water.

Phase	I	II	III	IV	V	VI
Step	Pre-rinsing	Cleaning	Neutralisation	Intermediate rinsing	Thermal disinfection	Drying
T (°C/°F)	< 25/77	55/131	20/68	70/158	94/201	90/194
t (min)	3	10	2	1	10	40
Water quality	D-W	P-W	P-W	P-W	P-W	-
Chemical specifications	-	Concentration, alkaline: pH = 10.9 < 5% anionic surfactants 1% solution pH = 10.5	Concentration, acidic: pH = 2.6; Base: Citric acid 1% solution pH = 3.0	-	-	-

D-W: drinking water, P-W: purified water demineralised

Checking, maintenance, testing

- Allow product to cool down to room temperature.
- Check product after each cleaning and disinfection for:
Cleanliness, functioning and damage, e.g. insulation, loose, bent, broken, cracked, worn and detached parts.
- Discard damaged products immediately. After the cleaning process, it is recommended to place the instrument inside suitable double sterile packaging. This will ensure reuse in the best way possible.

3.05 DISPOSAL

Adhere to national regulations when disposing of or recycling the product, its components or its packaging.

3.06 TROUBLESHOOTING

- The Verification Tool must be centred over the valve, otherwise the determined opening pressure may be incorrect.
- It must be ensured that the arrow markings on the Verification Tool correspond with the flow direction in the shunt system. The marking on the Verification Tool must face the proximal catheter (leading to the ventricle).
- When adjusting the adjustable Gravitational Unit, care must be taken that the opening pressure is changed by a maximum of 16 cmH₂O per adjustment process.
Example: The opening pressure for the adjustable Gravitational Unit on the proSA is to be changed from 6 to 36 cmH₂O. The correct method is an adjustment in two stages: Initial adjustment from 6 cmH₂O to 22 cmH₂O and then from 22 cmH₂O to 36 cmH₂O.

- ▶ The *proSA Verification Compass* reacts sensitively to external magnetic fields. In order to rule out unwanted interactions, the *proSA Adjustment Tool* should not be placed in the immediate vicinity of the *proSA Verification Compass* when determining the opening pressure. We recommend a minimum distance of 30 cm.
- ▶ Swelling of the skin may make adjustment difficult for a few days after the operation. If the valve setting cannot be checked conclusively using the proSA Tools, we recommend carrying out a check using an imaging method.
- ▶ The adjustable Miethke valve should not be implanted in an area that makes the detection or palpation of the valve difficult (e.g. underneath heavily scarred tissue).

3.07 TECHNICAL INFORMATION

3.07.01 TECHNICAL DATA

Manufacturer	Christoph Miethke GmbH & Co. KG
Product designation	<i>proSA Tools</i>
Specification of intended purpose	Treatment of hydrocephalus
Sterilisability	<i>proSA Checkmate</i> sterilisable <i>proSA Verification Tool</i> disinfectable <i>proSA Masterdisc</i> disinfectable <i>proSA Verification Compass</i> disinfectable <i>proSA Adjustment Tool</i> disinfectable <i>proSA Adjustment Disc</i> disinfectable
Storage	Store in a clean and dry place
Dimensions	<i>proSA Verification Tool</i> : ø 20 mm, L: 137 mm <i>proSA Masterdisc</i> : ø 39.5 mm, H: 11.6 mm <i>proSA Verification Compass</i> : ø 57.8 mm, H: 14.7 mm <i>proSA Adjustment Tool</i> : ø 20 mm, L: 137 mm <i>proSA Adjustment Disc</i> : ø 37.5 mm, H: 21.6 mm <i>proSA Checkmate</i> : ø 20 mm, H: 27.5 mm

3.07.02 MATERIALS THAT COME INTO CONTACT WITH BODILY TISSUE/FLUIDS

<i>proSA Verification Tool</i>	Brass (MS58), titanium (TiAl6V4)
<i>proSA Masterdisc</i>	Aluminium (AlSi1MgMn)
<i>proSA Verification Compass</i>	Aluminium (AlSi1MgMn)
<i>proSA Adjustment Tool</i>	Aluminium (AlMgSi0,5, AlMgSi1)
<i>proSA Adjustment Disc</i>	PPSU, PUR
<i>proSA Checkmate</i>	Titanium (TiAl6V4)

3.08 SYMBOLS USED FOR LABELLING

Symbol	Explanation
	EU conformity marking, xxxx indicates the number of the responsible Notified Body
	Medical device
	Manufacturer
	Date of manufacture
	Use by
	Production lot number, batch
	Article number
	UDI number (Unique Device Identifier)
	Steam-sterilised
	Do not resterilise
	Do not reuse
	Do not use if package is damaged
	Store in a dry place
	Temperature limitation
	Consult instructions for use / electronic instructions for use
	Caution, read the documentation supplied with the product
	Pyrogen-free
	Free of natural rubber latex, latex-free
	Indicates that in the USA, the product may only be issued to physicians

Symbol	Explanation
	MR unsafe
	Date
	Model number / European Medical Device Nomenclature Code

4.00 MEDICAL DEVICE CONSULTANTS

In compliance with regulatory requirements, Christoph Miethke GmbH & Co. KG has nominated medical device consultants as contacts for all product-related questions.

You can contact our medical device consultants at:

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany